### Compact Analytical Balances

# HIR-AZ/HIR-A

Series



HR-250AZ

HR-251AZ/HR-251A 62 g/252 g × 0.1 mg/1 mg

252 g x 0.1 mg

HR-250AZ/HR-250A

**HR-150AZ/HR-150A** 152 g x 0.1 mg

**HR-100AZ/HR-100A** 102 g x 0.1 mg



## Ultimate Value!





## Finally, Globally Proven Precision at a Fair Price

What do you expect most from a precision analytical balance? If you are like us, you will agree that accuracy and reliability come before anything else. The HR-AZ/HR-A series of compact analytical balances, behind their simple design, utilize proven A&D sensor technology that delivers superior performance to all users who wish to do superior work.

#### Optimizing value, the HR-AZ/HR-A series realizes the following three savings:

#### 1. Saving "Time"

#### Fast Stabilization in 2 Seconds

Our widely-recognized Compact Super Hybrid Sensor (C-SHS) technology now enables you to perform high-speed, frustration-free weighing even at a resolution of 0.1 mg. \*i Patent pending

#### One-Touch Automatic Calibration (HR-AZ models only)

The built-in calibration mass lets you ensure the accuracy of your balance with just one key press.

#### 2. Saving "Space"

#### Compact Footprint: 198 mm × 294 mm

With the C-SHS, the footprint of the balance is smaller than the ISO A4 (210 mm  $\times$  297 mm) paper size, maximizing opportunities for use of an analytical balance in even the smallest of spaces.

#### **Rotary Sliding Doors**

Unlike conventional analytical balances, the HR-AZ/HR-A series requires no extra space at the rear for accessing the weighing chamber, as the doors simply rotate behind the balance.

#### Removable Breeze Break

The breeze break is easily removed using a unique clip system, allowing fast, simple cleaning and use in confined spaces like gloveboxes and controlled environment cabinets.

**Rotary Sliding Doors** 



#### 3. Saving "Cost"

#### Lean Design for Affordability

The HR-AZ/HR-A design is all about "keeping it simple" without compromising quality or performance. Key features demanded by today's market have been included, resulting in a compact, precise, and reliable analytical balance for the value-conscious user.









































#### Superb Usability

#### Large Removable Breeze Break with Antistatic Coating

The breeze break is treated with an antistatic agent to ensure stable weighing. It also offers a large working space and can be detached for better portability and easy cleaning of the balance.



#### Clear, Reverse Backlit LCD

The contrast of the black and white display provides excellent readability even in a dimly-lit area.

#### Front-Mounted Spirit Level

The spirit level is easily visible when adjusting the leveling feet of the balance for accurate weighing.

#### **Great Flexibility and Compliance**

#### Twelve Weighing Units

You can easily choose from gram, milligram, ounce, troy ounce, metric carat, momme, pennyweight, grain, counting mode, percent mode, density mode, in addition to a user-programmable unit for conversion applications like GSM calculations in paper and fabric industries (one more unit from tael, tola, or Newton can be added upon request).

#### GLP/GMP/GCP/ISO Compliant Output

For compliance and documentation requirements, the HR-AZ/HR-A series can output the balance manufacturer, model, serial number, ID number, date and time, ispace for signature, and calibration records.

\*ii For HR-AZ models only. Please use the AD-8121B compact printer, which has a calendar function, to print the date and time with HR-A models.

#### Adjustable Response Characteristics and Hold Function

From three preset combinations of weighing speed and stability (FAST, MID, or SLOW), select the one that best suits your environmental conditions (drafts and vibrations). A hold function is also available for weighing animals.

#### Statistical Calculation Function

The balance displays and outputs statistical calculation data including number of data, sum, maximum, minimum, range (maximum–minimum), average, standard deviation, and coefficient of variation.

#### **Options and Accessories**

#### Quick USB Interface (FXi-02)

The easy-to-use USB interface requires no driver or software installation for transfer of weighing data to a PC.  $^{\rm iii}$ 

\*iii Data transfer from the balance to a PC only. Please use the RS-232C (standard interface) to serial-USB converter (AX-USB-9P, optional) for two-way communication or sending GLP or statistical calculation information via USB. Portable memory devices such as weighing environment logger (AD-1687) and data logger (AD-1688) are also available.

#### LAN-Ethernet Interface with WinCT-Plus Software (FXi-08)

You can send commands for controlling and acquiring data from multiple balances.



#### Rechargeable Battery Unit (FXi-09)

The Ni-MH rechargeable battery pack can be fitted into the balance for use without AC power (10 hours of charging for 14 hours of operation).

#### Density Determination Kit (AD-1654)

Using this kit, it is easy to calculate the density of a solid or liquid substance in density mode.

#### **Specifications**

		_								
Models		HR-250AZ	HR-251AZ	HR-150AZ	HR-100AZ	HR-250A	HR-251A	HR-150A	HR-100A	
Weighing capacity		252 g	62 g / 252 g iv	152 g	102 g	252 g	62 g / 252 g iv	152 g	102 g	
Minimum weighing value		0.1 mg	0.1 mg / 1 mg	0.1 mg		0.1 mg / 1 mg	0.1 mg			
Repeatability (standard deviation)		0 to 200 g: 0.1 mg	0.1 (0.5	0.1 mg		0 to 200 g: 0.1 mg	0.1 mg / 0.5 mg	0.1 mg		
		200 to 252 g: 0.2 mg	0.1 mg / 0.5 mg			200 to 252 g: 0.2 mg				
Linearity		±0.3 mg	±0.3 mg / ±1 mg	±0.2 mg ±0.3 mg		±0.3 mg / ±1 mg	±0.2 mg			
Stabilization time (when set to FAST under a good environment)		Approx. 2 seconds <sup>v</sup>								
Sensitivity drift		±2 ppm/°C (10 to 30 °C/50 to 86 °F)								
Internal calibration		Yes					No			
Clock and calendar function		Yes					No			
Operating environment		5 to 40 °C (41 to 104 °F), 85% RH or less (no condensation)								
Display refresh rate		5 times/second (10 times/second can be selected)								
Display mode Vi		gram, milligram, ounce, troy ounce, metric carat, momme, pennyweight, grain, counting mode, percent mode, density mode, and a u ser-programmable unit								
Counting mode	Minimum unit mass	0.1 mg	1 mg	0.1 mg			1 mg	0.1 mg		
	Number of samples	10, 25, 50 or 100 pieces								
Percent mode	Minimum 100% reference mass	10.0 mg	100 mg	10.0 mg			100 mg	10.0 mg		
	Minimum 100% display		0.01%, 0.1%, or 1% (depends on the reference mass stored)							
Interface		RS-232C								
Applicable calibration weight value		250 g, 200 g 100 g, 50 g	250 g, 200 g 100 g, 50 g, 20 g	150 g, 100 g, 50 g	100 g, 50 g	250 g, 200 g 100 g, 50 g	250 g, 200 g 100 g, 50 g, 20 g	150 g, 100 g, 50 g	100 g, 50 g	
Weighing pan size		Ø90 mm								
External dimensions		198 (W) × 294 (D) × 315 (H) mm								
Net weight		Approx. 3.9 kg Approx. 3.5 kg								
Power supply		AC adapter								
Power consumption		Approx. 11 VA (supplied to the AC adapter)								
¥ iv Smart range fur	ction: The minimum weighing val	uo will awitch to 1 may a	stamatically whom the d	ienlav valuo oveoode 63	a but ratures to 0.1 ma	humasina tha DE 7EE	O (tara) kau			

- 🛠 iv Smart range function: The minimum weighing value will switch to 1 mg automatically when the display value exceeds 62 g but returns to 0.1 mg by pressing the RE-ZERO (tare) key.
- \* v The factory setting is MID, whose stabilization time is approx. 3 seconds.
  vi One additional unit from tael (Singapore/HK jewelry/Taiwan/China), tola or Newton can be added upon request.

#### **Options**

FXi-02\* Quick USB interface with cable

FXi-08\* Ethernet interface with WinCT-Plus software

> Density determination kit Anti-vibration table for balances

FXi-09\* Built-in rechargeable battery

 $\ensuremath{\mbox{\textbf{\#}}}$  The FXi-02, FXi-08 and FXi-09 cannot be used at the same time

#### Accessories

AD-1654

AD-1671

AX-USB-9P


AD-1672	Tabletop breeze break
AD-1683	Static eliminator
AD-1684	Electrostatic field meter
AD-1687	Weighing environment logger
AD-1688	Weighing data logger
AD-1689	Tweezers for calibration weight
AD-8121B	Compact printer
AD-8526	Serial/Ethernet converter
AD-8527	Quick USB adapter
AD-8920A	Remote display
AD-8922A	Remote controller
AX-FXi-31	Main unit cover (5 pcs)


Serial/USB converter



(mm/inches)

**Dimensions** 









ACN 26 007 556 809 Head Office 32 Dew Street THEBARTON South Australia 5031 Telephone (08) 8301 8100 Facsimile (08) 8352 7409

Victorian Office 39 Bakehouse Road KENSINGTON Victoria 3031 Telephone (03) 9371 1555 Facsimile (03) 9372 1193

AD-1654

sales@andweighing.com.au www.andweighing.com.au

New South Wales Office Unit 4, 14 Abbott Road SEVEN HILLS New South Wales 2147 Telephone (02) 9674 5466 Facsimile (02) 9674 2544

Your Dealer