**Solvent Reclaimer
SR-305
Instruction Manual**

Preface

This instruction manual describes the operation of the Solvent Reclaimer, SR-305.

Be sure to read this manual before using the product to ensure proper and safe operation of the instrument.

Also safely store the manual so it is readily available whenever necessary.

Product specifications and appearance, as well as the contents of this manual are subject to change without notice.

■ Warranty and Responsibility

The product delivered to you is covered by HORIBA's warranty for a period of one (1) year.

If any malfunction or damage attributable to HORIBA's responsibility should occur during this period, necessary repairs or replacement of parts shall be made free of charge by HORIBA.

The warranty does not cover the following:

- Any malfunction attributable to improper operation
- Any malfunction attributable to repair or modification by any party not authorized by HORIBA
- Any malfunction attributable to the use in an improper environment
- Any malfunction attributable to violation of the instructions in this manual
- Any malfunction attributable to operations in the manner not specified in this manual
- Any malfunction attributable to natural disasters, or accidents or mishaps not involving HORIBA
- Any deterioration in appearance attributable to corrosion, rust, and so on.
- Consumables and replacement of consumables
- Products of other companies

HORIBA shall not be liable for any damages resulting from any malfunctions of this product, any erasure of data, or any other uses of this product.


■ Trademarks

Generally, company names and brand names are either registered trademarks or trademarks of the respective companies.

Safety Notice

■ The Safety Messages Used in Manual and Their Meanings

The safety alerts that are used in this manual consist of a warning symbol, and symbol word. You must obey these safety alerts.

● Signal words and their meanings

- CAUTION:
Indicates latent danger. If the danger is not avoided personal injury may result. This alert is also used to indicate unsafe behavior.

■ General Points to Note Regarding the SR-305

- The SR-305 is the exclusive solvent reclaimer for S-316 and/or H-997. Do not use the SR-305 for other solvents.
- The solvent (H-997) is listed on Montreal Protocol Annex C Group I. For its consumption and re-export, strict obedience to your local laws is necessary.
- Carefully read the check and maintenance section in this manual before attempting maintenance work.
- To avoid risk of electric shock, never remove any covers from the SR-305 unless specified in this manual.

■ If You Do Not Follow the Procedures and Cautions

- HORIBA assumes no responsibility for matters which occur due to a failure to follow the instructions and WARNINGS described in this manual.

■ Handling Cautions

- Avoid operating or storing the SR-305 in the following locations or situations:
 - Where the humidity is above 80%.
 - Where the temperature is below 0°C or above 40°C.
 - Where it will be exposed to direct sunlight.
 - Dusty locations.
 - Poor ventilation.
 - Presence of corrosive gases.
- Use a protective mask and polyethylene gloves under good ventilation when handling solvents.
- Background readings may differ slightly depending on the solvents manufacturing lot. When using solvents from different manufacturing lots or reclaimed solvent, match the backgrounds by mixing the amount of solvent that you plan to use in a glass container before using it.
- The SR-305 is the exclusive solvent reclaimer for S-316 and/or H-997.
Do not use the SR-305 for other solvents.

■ Structure of This Manual

This manual is composed of the 5 chapters listed below.

Chapter 1. Checking Contents

Shows a set of the SR-305.

Chapter 2. Parts Description

Gives an outline of the SR-305.

Chapter 3. Preparation

Describes the preparation before reclaiming for the SR-305.

Chapter 4. Operating Procedures

Describes the reclaiming procedures for the SR-305.

Chapter 5. Regular Maintenance

Describes the daily maintenance procedures for the SR-305.


Contents

1	Checking Contents	1
2	Parts Description	2
3	Setting Up the SR-305	4
4	Operation Procedure	6
4.1	Reclamation of S-316	6
4.2	Reclamation of H-997	7
5	Daily Maintenance	8
5.1	Replacing active carbon and alumina	8
5.2	Replacing the water separation filter	8
5.3	Spare parts	8

1 Checking Contents

The package includes the following main unit and accessories. Check that none of the items listed below is missing or damaged.


Main unit - 1 set

Aluminium oxide active acidic -1 bottle (90 g)


Water separating filter (120 mm dia.) - 1 sheet

Activated carbon filter (100 mm dia.) - 1 sheet

Activated carbon - 1 can (500 g)


O-ring - 1 pc



Instruction manual - 1 copy

2 Parts Description

● Dimensional outline


[Unit: mm]


● Mass: Approx. 7 kg

3 Setting Up the SR-305


1. Fit the water separation filter (the larger filter) to the water separation tank.


2. Pull the latch to remove the active carbon tank and the alumina tank.


3. Fill the alumina tank with all of alumina (90 g).

-
4. Fit the active carbon filter (the smaller filter) to the active carbon tank, and fill the tank with all of active carbon (500 g). Put the tank back in place and lock it with the latch.

5. Put the empty bottle in place, which is used to reserve the reclaimed solvent.

4 Operation Procedure

4.1 Reclamation of S-316

1. Use the calibrated OCMA-300 to measure the concentration of the solvent being reclaimed.
2. Open the lid of the SR-305 and put the solvent. A small amount of water that may exist in the water separation vessel (up to about 1 cm from the bottom) is negligible, but frequently drain water from the vessel.

CHEMICAL HAZARD (S-316)

S-316 can be absorbed through the skin and a large amount can cause poisoning. Wear gloves when handling it, and take precautions to ensure that you never ingest it accidentally.

Note

One cycle of reclamation fills up the empty bottle (about 500 mL).

- When active carbon is new, the solvent is absorbed by active carbon, the recovery efficiency of solvent is about 35% in the first cycle. It should be about 90% in each of the subsequent cycles.
 - When active carbon is new, the active carbon vessel generates heat with a reaction during the reclamation, but this is not abnormal.
 - Close the lid and the door during the reclamation.
-

3. The reclamation occurs at a rate of about 4 mL/min to 6 mL/min, which means that it takes about 120 minutes to 180 minutes to fill up each bottle.
4. Use the calibrated OCMA to measure the concentration of the solvent reclaimed.
5. If that concentration is no lower than 5 mg/L, reclaim the solvent again to prevent the concentration of oil contents in the solvent from increasing due to the insufficiency of reclaiming the solvent.

Note

- The solvent may not be completely reclaimed in one cycle depending on its oil type such as OCB and its concentration.
 - Any solvent containing surface active agent (an emulsifying agent layer) deteriorates the function of the water separation filter. In this case, remove the surface active agent (emulsifying agent layer) and then reclaim the solvent.
 - If water is reserved in the water separation vessel exceeding about 1 cm from the bottom, the water may leak from the water separation filter. Therefore, frequently drain water from the vessel. To do this, open the lid and remove the tank to drain water.
-

Reference

For a span calibration solution prepared with OCB, the concentration of oil contents is more difficult to decrease. One reclamation cycle allows for decreasing that concentration to about 70% of the previous value. The concentration of oil contents can be decreased to the desired level by repeating the reclamation cycle. Two or three reclamation cycles allows for decreasing the concentration to below 5 mg/L. If the initial concentration of oil contents is no lower than 20 mg/L, more than two reclamation cycles are required. For this reason, it is recommended that you reclaim such a solvent separately from the solvents which contain other kinds of oil.

4.2 Reclamation of H-997

● Preparation

- Settle the SR-305 in a place with good ventilation such as a draft chamber for reclaiming H-997.
- Prepare over 600 mL of used solvent to be reclaimed.
- Remove water from the solvent with pipette as much as possible.
- Measure the oil concentration in the solvent with calibrated OCMA-3X5 in advance.

● Operation

CAUTION

CHEMICAL HAZARD (H-997)

Ingestion and Inhalation of a large amount of H-997 can cause poisoning. Please take precautions to ensure that you have not ingested or inhaled it accidentally.

1. Carefully put 300 mL solvent to be reclaimed in the water separation tank on the top of the SR-305.
2. In case the solvent passes the new active carbon, the carbon generates heat. Leave the reclaimer for over 30 minutes until the carbon cools down to the room temperature.

Note

Temperature of the active carbon rises to approximately 70°C, but this usually happens. The solvent is hardly recovered, since the solvent is adsorbed and vaporizes during passing through the heated active carbon. Close the lid and the door during the reclamation.

3. Put the solvent (300 mL to 500 mL) in the water separation tank again.
4. Measure the oil concentration of recovered solvent with calibrated OCMA-3X5, and reuse that for further analysis if the concentration is below 5 mg/L. Start the reclamation from step 3. again, if the measured concentration is over 5 mg/L.
5. For reclaiming the solvent by the active carbon which the solvent passed through once, start the reclamation from step 3.

Reference

- In most cases, the active carbon does not heat after the carbon experienced the contact to the solvent once, and volumetric recovery rises to 60%. One time reclamation enables that over 98% (mass) of oil dissolved in solvent will be removed. For example, when the used solvent (500 mL) which contains 25 mg/L of oil concentration is reclaimed, volumetric recovery of the solvent will be over 300 mL and its oil concentration will drop to less than 1 mg/L.
- It takes 30 minutes to 45 minutes to reclaim 500 mL solvent.
- The solvent may not be completely reclaimed in one reclamation depending on its oil type such as OCB and its concentration.
- Any solvent containing surfactant (emulsified layer) spoils the function of the water separation filter. Remove the emulsified layer from the solvent in advance.
- Drain trapped water in the water separation tank frequently, because the water may leak through the water separation filter if water level exceeds 1 cm from the filter surface. (The water separation tank can be taken out from the top of SR-305.)
- After the reclamation has finished, it is recommended that the active carbon is transferred in a sealed vessel and stored in a dark and cool room to prevent the solvent vaporized from the carbon.

5 Daily Maintenance

5.1 Replacing active carbon and alumina

- **For reclamation of S-316**

The SR-305 has an adsorption capacity of 1000 mg. Use the following equation to determine the frequencies of replacing active carbon and alumina from the concentration of oil contents in your sample and total amount of reclamation.

$$1000 \text{ (mg)} / \text{concentration of oil contents in the solvent being reclaimed (mg/L)} \\ = \text{reclaimed amount (L)}$$

e.g.: A solvent whose concentration of oil contents is 20 mg/L can be reclaimed up to an amount of 50 L.

Replace active carbon alumina using the same procedure as described in " 3 Setting Up the SR-305 " (page 4) The active carbon filter must be fit after it is well dried with air. If the active carbon filter is broken, replace it.

- **For reclamation of H-997**

The SR-305 has an adsorption capacity of 1400 mg. Use the following equation to determine the frequencies of replacing active carbon and alumina from the concentration of oil contents in your sample and total amount of reclamation.

$$1400 \text{ (mg)} / \text{concentration of oil contents in the solvent being reclaimed (mg/L)} \\ = \text{reclaimed amount (L)}$$

e.g.: A solvent whose concentration of oil contents is 20 mg/L can be reclaimed up to an amount of 70 L.

Replace active carbon alumina following the same procedure as described in " 3 Setting Up the SR-305 " (page 4) The active carbon filter must be fit after it is well dried with air. If the active carbon filter is broken, replace it.

5.2 Replacing the water separation filter

When the water separation filter become dirty or reclamation rate becomes less than 6 mL/min, replace the water separation filter. To do this, follow the same procedure as described in " 3 Setting Up the SR-305 " (page 4).

5.3 Spare parts

Part Name	Part Number	Specifications
O-Ring	9105000100	Contains 2 pcs. of type P110
Water Separation Filter	9039001800	Contains 2 pcs. of type 120 mm dia.
Active Carbon Filter	9039001700	Contains 2 pcs. of type 100 mm dia.
Active Carbon	9039000500	Contains 500 g of type BAC-MP
Active Alumina	9039001900	Contains 90 g

These supplies may be ordered by part name and part number from your dealer.

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© 2005 HORIBA, Ltd.