**Guardian™ 2000**  
**Hotplate-Stirrer, e-G21**  
**Instruction Manual**

## Table of Contents

1.	INTRODUCTION .....	2
1.1.	Safety Information.....	2
1.2.	Intended Use .....	3
1.3.	Package Contents .....	3
1.4.	Installation.....	4
2.	OVERVIEW .....	5
2.1.	Dimensions .....	5
2.2.	Control Panel .....	8
3.	OPERATION .....	9
3.1.	Getting Ready .....	9
3.2.	Standby Mode.....	9
3.3.	Controlling the Stirrer .....	9
3.4.	Controlling the Top Plate Heater.....	10
4.	MAINTENANCE .....	11
4.1	Cleaning.....	11
4.2	Replacing Power Fuse.....	11
4.3	Troubleshooting .....	13
4.4	Service Information.....	13
5.	TECHNICAL DATA.....	14
6.	COMPLIANCE.....	14

## 1. INTRODUCTION

This manual contains installation, operation and maintenance instructions for the OHAUS Guardian™ 2000 hotplate-stirrer. Please read the manual completely before using.

### 1.1. Safety Information

#### Definition of Signal Warnings and Symbols

Safety notes are marked with signal words and warning symbols. These show safety issues and warnings. Ignoring the safety notes may lead to personal injury, damage to the instrument, malfunctions and false results.

WARNING	For a hazardous situation with medium risk, possibly resulting in severe injuries or death if not avoided.
CAUTION	For a hazardous situation with low risk, resulting in damage to the device or the property or in loss of data, or minor or medium injuries if not avoided.
ATTENTION	For important information about the product. May lead to equipment damage if not avoided.
NOTE	For useful information about the product.

#### Warning Symbols

General hazard

Explosion hazard

Caution, hot surface

Protective conductor terminal

Alternating current

Electrical shock hazard

#### Safety Precautions

**WARNING!** The protection provided by the unit may be impaired if used with accessories not provided or recommended by the manufacturer, or used in a manner not specified by the manufacturer.

- Always operate unit on a level surface for best performance and maximum safety.
- DO NOT lift unit by the top plate.
- To avoid electrical shock, completely cut off power to the unit by disconnecting the power cord from the unit or unplugging from the wall outlet.
- Disconnect unit from the power supply prior to maintenance and servicing.
- Spills should be removed promptly, after the unit has cooled down.
- Alkalis spills, hydrofluoric acid or phosphoric acid spills may damage the unit and lead to thermal failure.
- DO NOT immerse the unit for cleaning.
- DO NOT operate the unit at high temperatures without a vessel/sample on the top plate.
- DO NOT operate the unit if it shows signs of electrical or mechanical damage.
- Protective earthing of the equipment is achieved via connection of the provided power cord to a compatible grounded power outlet.

**WARNING:** unit is not explosion proof. Use caution when unit is on or when heating volatile materials.

**WARNING!** DO NOT use the unit in explosive atmospheres or with materials that could cause a hazardous environment from processing. Keep in mind the material flash point relative to the target temperature that has been set.

**CAUTION!** The top plate can reach 500°C, DO NOT touch the heated surface. Use caution at all times. Keep the unit away from explosive vapors and clear of papers, drapery and other flammable materials. Keep the power cord away from the heater plate.

**CAUTION! Beware of the following risks when heating.**

- Flammable materials
- Low boiling point combustible substances
- Glass breakage as a result of mechanical shaking power
- Incorrect container size
- Too much medium
- Unsafe condition of container

Earth Ground - Protective Conductor Terminal

Alternating Current

## 1.2. Intended Use

This instrument is intended for use in laboratories, pharmacies, schools, businesses and light industry. It must only be used for processing materials as described in these operating instructions. Any other type of use and operation beyond the limits of technical specifications, without written consent from OHAUS, is considered as not intended. This instrument complies with current industry standards and the recognized safety regulations; however, it can constitute a hazard in use. If the instrument is not used according to these operating instructions, the intended protection provided by the instrument may be impaired.

## 1.3. Package Contents

- Hotplate Stirrer
- Power Cord
- Stir Bar

## 1.4. Installation

Upon receiving the Ohaus Hotplate-Stirrer, check to ensure that no damage has occurred during shipment. It is important that any damage that occurred in transport is detected at the time of unpacking. If you do find such damage, the carrier must be notified immediately.

After unpacking, place the Hotplate-Stirrer on a level bench or table, away from explosive vapors. Ensure that the surface on which the unit is placed will withstand typical heat produced by the unit and place the unit a minimum of six (6) inches from vertical surfaces. Do not position the equipment such that it is difficult to disconnect the power cord during use. Always place the unit on a sturdy work surface.

The Hotplate-Stirrer is supplied with a 3 conductor, grounded power cord that should be plugged into a matching standard grounded outlet. If the cord supplied does not meet your needs, please use an approved power cord that has ratings equal or exceeding those of the originally provided cord and that complies with the local/national regulations of the country in which the equipment is to be used. Replacement of the plug must be made by a qualified electrician.

## 2. OVERVIEW

### 2.1. Dimensions

#### Round Top Hotplate-Stirrer

Overall Dimensions (L x W x H)	264 x 163 x 101 mm (10.4 x 6.4 x 4.0")
Top Plate Dimensions	Ø 13.5 cm (5.3")
Top Plate Material	Ceramic coated stainless steel
Electrical (50 / 60 Hz) - Hotplate Stirrer - Hotplate	100 - 120 V ~ 5 A / 220 - 240 V ~ 4 A
Electrical (50 / 60 Hz) - Stirrer	100 - 120 V ~ 2 A / 220 - 240 V ~ 1 A
Fuses - Hotplate Stirrer - Hotplate	10 A, 5 x 20 mm, 250 VAC / 6.3 A, 5 x 20 mm, 250 VAC
Fuses - Stirrer	6.3 A, 5x20 mm, 250 VAC
Temperature Range	50° to 380°C
Temperature Stability of Top Plate <sup>+</sup>	± 3% (> 100° C), ± 3° C (≤ 100° C)
Stir Capacity	15 L
Speed Range <sup>**</sup>	200 to 1600 rpm
Weight Capacity	Up to 15 kg (33 lbs)
Shipping Weight - Hotplate Stirrer - Hotplate - Stirrer	3.1 kg 2.4 kg 2.9 kg

#### Note:

- <sup>+</sup> The parameter in the specification table are applicable to 2" (5 cm) diameter center of the top plate.
- <sup>\*\*</sup> Minimal speed varies depending on the stir bar, vessel and fluid conditions.

## 4x4 Hotplate-Stirrer

Overall Dimensions (L x W x H)	264 x 163 x 110 mm (10.4 x 6.4 x 4.3")
Top Plate Dimensions	10.2 x 10.2 cm (4 x 4")
Top Plate Material - Hotplate Stirrer - Hotplate	Ceramic
Top Plate Material - Stirrer	Resin
Electrical (50 / 60 Hz) - Hotplate Stirrer - Hotplate	100 – 120 V ~ 5 A / 220 – 240 V ~ 4 A
Electrical (50 / 60 Hz) - Stirrer	100 – 120 V ~ 2 A / 220 – 240 V ~ 1 A
Fuses - Hotplate Stirrer - Hotplate	10 A, 5 x 20 mm, 250 VAC / 6.3 A, 5 x 20 mm, 250 VAC
Fuses - Stirrer	6.3 A, 5x20 mm, 250 VAC
Temperature Range	70° to 500°C
Temperature Stability of Top Plate <sup>+</sup>	± 3% (> 100° C), ± 3° C (≤ 100° C)
Stir Capacity	15 L
Speed Range <sup>++</sup>	200 to 1600 rpm
Weight Capacity	Up to 15 kg (33 lbs)
Shipping Weight - Hotplate Stirrer - Hotplate - Stirrer	3.1 kg 2.4 kg 2.6 kg

## Note:

- <sup>+</sup> The parameter in the specification table are applicable to 2" (5 cm) diameter center of the top plate.
- <sup>++</sup> Minimal speed varies depending on the stir bar, vessel and fluid conditions.

## 7x7 Hotplate-Stirrer


Overall Dimensions (L x W x H)	318 x 223 x 112 mm (12.5 x 8.8 x 4.4")
Top Plate Dimensions	17.8 x 17.8 cm (7 x 7")
Top Plate Material - Hotplate Stirrer - Hotplate	Ceramic
Top Plate Material - Stirrer	Resin
Electrical (50 / 60 Hz) - Hotplate Stirrer - Hotplate	100 - 120 V ~ 10 A / 220 - 240 V ~ 7 A
Electrical (50 / 60 Hz) - Stirrer	100 - 120 V ~ 2 A / 220 - 240 V ~ 1 A
Fuses - Hotplate Stirrer - Hotplate	15 A, 5 x 20 mm, 250 VAC / 10 A, 5 x 20 mm, 250 VAC
Fuses - Stirrer	6.3 A, 5x20 mm, 250 VAC
Temperature Range	70° to 500°C
Temperature Stability of Top Plate <sup>+</sup>	± 3% (> 100° C), ± 3° C (≤ 100° C)
Stir Capacity	15 L
Speed Range <sup>++</sup>	200 to 1600 rpm
Weight Capacity	Up to 15 kg (33 lbs)
Shipping Weight - Hotplate Stirrer - Hotplate - Stirrer	4.6 kg 3.9 kg 3.5 kg

## Note:

- <sup>+</sup> The parameter in the specification table are applicable to 2" (5 cm) diameter center of the top plate.
- <sup>++</sup> Minimal speed varies depending on the stir bar, vessel and fluid conditions.

## 2.2. Control Panel

1. **Left Knob:** Controls temperature.
2. **Standby Switch**
3. **Heater Indicator:** Illuminates when heater is running.
4. **Hot Top Caution Indicator:** Illuminates when the top plate surface is hot. Do not touch!
5. **Speed Indicator:** Illuminates when stirrer is running.
6. **Right Knob:** Controls speed.
7. **Power Entry Module (PEM)**
8. **Fuse**
9. **Accessory Rod Mounting Hole**

**Note:** The top plates are of different dimensions and materials. Please refer to the actual product for details.

### 3. OPERATION

#### 3.1. Getting Ready

To get ready:

1. Plug the female end of the provided power cord into **Power Entry Module (PEM)** on the rear side of the unit.
2. Plug the male end of the power cord into a matching standard grounded outlet.

#### 3.2. Standby Mode

The standby switch on the front side of the unit controls standby mode.

1. When the unit is switched off:
  - a) All heating and stirring function will turn off.
  - b) If the top plate is hot, the  on the control panel will remain illuminated.
2. When the unit is switched on:
  - a) The unit is ready for normal use.
  - b) If the heating and stirring knobs are turned off, heating and stirring functions will remain off.
  - c) If the heating or stirring knobs are turned on, it starts heating or stirring immediately.

#### 3.3. Controlling the Stirrer

1. Rotate the right knob to control the speed setting.
  - a) Clockwise rotation will increase the speed setting.

- b) Counterclockwise rotation will decrease the speed setting.

- To turn on the stirrer, rotate the knob clockwise from the OFF position. The  illuminates to indicate that the stirrer is on.
- To turn off the stirrer, rotate the knob counterclockwise and return to the OFF position. The  light turns off to indicate that the stirrer is off.

### STIRRING OPERATING TIPS

The stirrer increases speed at a steady rate until the setpoint is reached. If the stirrer is not reaching its setpoint: 1) the stir bar may be too large, 2) the liquid may be too viscous, 3) the setpoint speed may need to be reduced. Additionally, the magnetic strength of stir bars reduce over time and may need to be replaced. Apart from that, the bottom surface of the vessel can affect the stirring performance of the magnetic stir bar. Adjust the vessel being used or change the stir bar to optimize mixing in the sample.

When heating and stirring a reaction vessel within an oil bath or similar set-up, the stirring function will stir up to approximately one inch (2.5 cm) from the top plate. The stirring speed will vary according to liquid viscosity, spin bar length, and distance from top plate. Adjust one or all of these to achieve the desired stirring speed. For example: the closer the reaction vessel is to the top plate, the stronger the magnetic connection between the unit and the stir bar.

### 3.4. Controlling the Top Plate Heater

- Rotate the left knob to control the heat setting.
  - Clockwise rotation will increase the heat setting.

- Counterclockwise rotation will decrease the heat setting.

- To turn on the heater, rotate the knob clockwise from the OFF position. The  illuminates to indicate that the heater is running.
- To turn off the heater, rotate the knob counterclockwise and return to the OFF position. The  light turns off to indicate that the heater is off.

**CAUTION: THIS DOES NOT MEAN THAT THE TOP PLATE IS SAFE TO TOUCH.**

Once the heater is turned off and the heater's temperature cools, the  on the control panel will disappear.

**Typical Time to Boil Water**

The chart below is an example of an approximate time to boil for the specified amount of water in a specific vessel. These values are only approximate and can vary from unit to unit. Values are based on 23°C water in an ambient environment of 23°C.

Unit Size	Heater Temp. Limit	Volume of Water	Typical Time to Boil
Round Top	380°C	500mL in 1L beaker	≈ 25 min
7×7	500°C	500mL in 1L beaker	≈ 18 min
4×4	500°C	500mL in 1L beaker	≈ 30 min

**4. MAINTENANCE****4.1 Cleaning**

**WARNING:** Electric Shock Hazard. Disconnect the equipment from the power supply before cleaning.

Make sure that no liquid enters the interior of the instrument.


**Attention:** Do not use solvents, harsh chemicals, ammonia or abrasive cleaning agents.

The housing may be cleaned with a cloth dampened with a mild detergent if necessary.

**4.2 Replacing Power Fuse**


The use of a fuse of a different type or with a different value, or bridging or shunting the fuse is not allowed and can possibly cause a hazard to your safety and lead to instrument damage!


**WARNING: Electric Shock Hazard.** Disconnect the equipment from the power supply before replacing the fuse.

To replace the fuse:

1. Put the tip of a screwdriver or a test pen in the circled position below, and then pull the fuse holder out.


2. Take out the old fuse and then push the new fuse into the slot. Choose the correct fuse to replace according to the fuse specification listed in the **Dimensions** Section.

**Note:** the unit comes with a spare fuse in the fuse holder (in the circled position above).

**Attention:** If the fuse is good and power is available at the outlet, the cord or unit may be defective. Try a new cord. If this does not work, the unit should be sent back for servicing.

3. Push the fuse holder back in to finish the replacement.

### 4.3 Troubleshooting

The following table lists common problems and possible causes and remedies. If the problem persists, contact OHAUS or your authorized dealer.

Table 4-1 Troubleshooting

<b>Error</b>	<b>Cause of Error</b>	<b>How to Fix</b>
Unit fails to power on	Missing or blown fuse	Add or replace fuse as necessary.
Unit fails to heat	Heater over temperature or Heater damage	Unplug and replug the power cord for reset, if the failure consistent, please contact Ohaus.
Unit fails to stir	Motor overload or Motor damage	Smaller vessel or reduce volume in vessel and wait for auto reset, if the failure consistent, please contact Ohaus.

### 4.4 Service Information

If the troubleshooting section does not resolve your problem, contact an authorized OHAUS Service Agent. For Service assistance in the United States, call toll-free 1-800-526-0659 between 8:00 AM and 5:00 PM Eastern Standard Time. An OHAUS Product Service Specialist will be available to assist you. Outside the USA, please visit our website [www.ohaus.com](http://www.ohaus.com) to locate the OHAUS office nearest you.

## 5. TECHNICAL DATA

Equipment Ratings:

Indoor use only

Altitude: 0 to 2000m

Operating temperature: 5°C to 40°C

Storage temperature: -20°C to 65°C

Operating humidity: 20 to 80% relative humidity, non-condensing

Storage humidity: 20 to 80% relative humidity, non-condensing

Voltage fluctuations: Mains supply voltage fluctuations up to  $\pm 10\%$  of the nominal voltage.

Overvoltage category II

(Installation category):

Pollution degree: 2

## 6. COMPLIANCE

Compliance to the following standards is indicated by the corresponding mark on the product.

Mark	Standard
	This product complies with the applicable harmonized standards of EU Directives 2011/65/EU (RoHS), 2014/30/EU (EMC) and 2014/35/EU (LVD). The EU Declaration of Conformity is available online at <a href="http://www.ohaus.com/ce">www.ohaus.com/ce</a> .
	This product complies with the applicable statutory standards of the Restriction of the Use of Certain Hazardous Substances in Electrical and Electronic Equipment Regulations 2012, UK Electromagnetic Compatibility Regulations 2016 and Electrical Equipment (Safety) Regulations 2016. The UK Declaration of Conformity is available online at <a href="http://www.ohaus.com/uk-declarations">www.ohaus.com/uk-declarations</a> .
	This product complies with the EU Directive 2012/19/EU (WEEE). Please dispose of this product in accordance with local regulations at the collecting point specified for electrical and electronic equipment. For disposal instructions in Europe, refer to <a href="http://www.ohaus.com/weee">www.ohaus.com/weee</a> .
	EN 61326-1
	CAN/CSA-C22.2 No. 61010-1, CAN/CSA-C22.2 No. 61010-2-010, CAN/CSA-C22.2 No. 61010-2-051
	UL 61010-1, UL 61010-2-010, UL 61010-2-051

### ISED Canada Compliance Statement:

CAN ICES-003(A) / NMB-003(A)

### ISO 9001 Registration

The management system governing the production of this product is ISO 9001 certified.

**FCC Supplier Declaration of Conformity**

Unintentional Radiator per 47CFR Part B  
Trade Name: OHAUS CORPORATION  
Model or Family identification: e-G21xxxxx

**Party issuing Supplier's Declaration of Conformity:**

Ohaus Instruments (Changzhou) Co., Ltd.  
2F, 22 Block, 538 West Hehai Road, Xinbei District, Changzhou  
Jiangsu 213022  
China  
Phone: +86 519 85287270

**Responsible Party – U.S. Contact Information:**

Ohaus Corporation  
7 Campus Drive, Suite 310  
Parsippany, NJ 07054  
United States  
Phone: +1 973 377 9000  
Web: [www.ohaus.com](http://www.ohaus.com)

**FCC Compliance Statement:**

Note: This equipment has been tested and found to comply with the limits for a Class A digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference when the equipment is operated in a commercial environment. This equipment generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with the instruction manual, may cause harmful interference to radio communications. Operation of this equipment in a residential area is likely to cause harmful interference in which case the user will be required to correct the interference at his own expense.

Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

**LIMITED WARRANTY**

OHAUS products are warranted against defects in materials and workmanship from the date of delivery through the duration of the warranty period. During the warranty period OHAUS will repair, or, at its option, replace any component(s) that proves to be defective at no charge, provided that the product is returned, freight prepaid, to OHAUS.

This warranty does not apply if the product has been damaged by accident or misuse, exposed to radioactive or corrosive materials, has foreign material penetrating to the inside of the product, or as a result of service or modification by other than OHAUS. In lieu of a properly returned warranty registration card, the warranty period shall begin on the date of shipment to the authorized dealer. No other express or implied warranty is given by OHAUS Corporation. OHAUS Corporation shall not be liable for any consequential damages.

As warranty legislation differs from state to state and country to country, please contact OHAUS or your local OHAUS dealer for further details.

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P/N 30692776 A © 2022 Ohaus Corporation, all rights reserved / todos los derechos reservados / tous droits réservés / alle Rechte vorbehalten / tutti i diritti riservati