Instruction Manual

pH/ION/COND METER F-74BW

Preface

This manual describes the operation of the F-74BW pH/ION/COND Meter. Be sure to read this manual before using the product to ensure proper and safe operation of the instrument. Also safely store the manual so it is readily available whenever necessary.

Product specifications and appearance, as well as the contents of this manual are subject to change without notice.

Warranty and Responsibility

HORIBA, Ltd. warrants that the Product shall be free from defects in material and workmanship and agrees to repair or replace free of charge, at option of HORIBA, Ltd., any malfunctioned or damaged Product attributable to responsibility of HORIBA, Ltd. for a period of one (1) year from the delivery unless otherwise agreed with a written agreement. In any one of the following cases, none of the warranties set forth herein shall be extended;

- Any malfunction or damage attributable to improper operation
- Any malfunction attributable to repair or modification by any person not authorized by HORIBA, Ltd.
- Any malfunction or damage attributable to the use in an environment not specified in this manual
- Any malfunction or damage attributable to violation of the instructions in this manual or operations in the manner not specified in this manual
- Any malfunction or damage attributable to any cause or causes beyond the reasonable control of HORIBA, Ltd. such as natural disasters

• Any deterioration in appearance attributable to corrosion, rust, and so on Replacement of consumables

HORIBA, LTD. SHALL NOT BE LIABLE FOR ANY DAMAGES RESULTING FROM ANY MALFUNCTIONS OF THE PRODUCT, ANY ERASURE OF DATA, OR ANY OTHER USES OF THE PRODUCT.

Trademarks

Generally, company names and brand names are either registered trademarks or trademarks of the respective companies.

Unauthorized Reprinting or Copying of This Operation Manual

No unauthorized reprinting or copying of all or part of this operation manual is allowed. The utmost care has been used in the preparation of this operation manual.

Copyright © 2011 by HORIBA, Ltd. All rights reserved.

Conformable Directive

This equipment conforms to the following directives and standards:

CE

Directives: The EMC Directive 2004/108/EC The Low Voltage Directive 2006/95/EC
Standards: [the EMC Directive] EN61326-1:2006 Class B, Basic requirements

[the Low Voltage Directive] EN61010-1:2010(Ed.3.0)

Installation Environment

This product is designed for the following environment.

- Overvoltage category II
- Measurement category I

WARNING: Do not use the equipment for measurements within measurement categories II, III and IV.


Information on Disposal of Electrical and Electronic Equipment and Disposal of Batteries and Accumulators

The crossed out wheeled bin symbol with underbar shown on the product or accompanying documents indicates the product requires appropriate treatment, collection and recycle for waste electrical and electronic equipment (WEEE) under the Directive 2002/96/EC, and/or waste batteries and accumulators under the Directive 2006/66/EC in the European Union.

The symbol might be put with one of the chemical symbols below. In this case, it satisfies the requirements of the Directive 2006/66/EC for the object chemical.

This product should not be disposed of as unsorted household waste. Your correct disposal of WEEE, waste batteries and accumulators will contribute to reducing wasteful consumption of natural resources, and protecting human health and the environment from potential negative effects caused by hazardous substance in products.

Contact your supplier for information on applicable disposal methods.



FCC Rules

Any changes or modifications not expressly approved by the party responsible for compliance shall void the user's authority to operate the equipment.

WARNING

This equipment has been tested and found to comply with the limits for a Class A digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference when the equipment is operated in a commercial environment. This equipment generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with the instruction manual, may cause harmful interference to radio communications.

Operation of this equipment in a residential area is likely to cause harmful interference in which case the user will be required to correct the interference at his own expense.

Hazard Classification and Warning Symbols

Warning messages are described in the following manner. Read the messages and follow the instructions carefully.

Hazard classification

This indicates an imminently hazardous situation which, if not ▲ DANGER avoided, will result in death or serious injury. This signal word is to be limited to the most extreme situations.



A CAUTION

MARNING This indicates a potentially hazardous situation which, if not avoided, could result in death or serious injury.

This indicates a potentially hazardous situation which, if not avoided, may result in minor or moderate injury. It may also be used to alert against unsafe practices. Without safety alert indication of hazardous situation which, if not avoided, could result in property damage.

Warning symbols



Description of what should be done, or what should be followed



Description of what should never be done, or what is prohibited

Safety Precautions

This section provides precautions to enable you to use the product safely and correctly and to prevent injury and damage. The terms of DANGER, WARNING, and CAUTION indicate the degree of imminency and hazardous situation. Read the precautions carefully as it contains important safety messages.

0	Do not use an unspecified AC adapter. Otherwise, it may heat up or be ignited resulting in a fire or an accident.
0	Do not disassemble or modify the meter. Otherwise, it may heat up or be ignited resulting in a fire or an accident.
	CAUTION
9	Harmful chemicals Some ion electrodes are used with hazardous standard solutions. Handle them with care. If the internal solution comes in contact with the skin, wash it off immediately. If it gets into eyes, flush with plenty of water and then consult a doctor.
0	Harmful chemicals The internal solution of an electrodes is highly concentrated potassium chloride (3.33 mol/L KCl). If the internal solution comes in contact with the skin, wash it off immediately. If it gets into eyes, flush with plenty of water and then consult a doctor.
0	Broken glass Broken glass may cause injury. The outer tube and tip of an electrode are made of glass. Handle them with care.
0	Do not use the cable of serial communication, USB, or AC adapter under wet or humid conditions. Otherwise, it may cause an fire, electric shock, or breakage.

Product Handling Information

Operational Precautions

- Do not drop, crash, or give any physical impact on the instrument.
- Do not immerse the instrument into alcohol, organic solvent, strong acid, strong alkaline, or the like. The instrument body contains ABS resin, acrylic resin, and some rubber parts.
- If the instrument is dropped into water or gets wet, wipe it using soft cloth. Do not heat to dry it with a hair-dryer (or the like).
- Use fingers to press the operation keys or the touch panel. Do not use a hard object like a metal stick or rod.
- Be careful not to let water into the instruction inside. The instrument is not water-proof.
- To disconnect an electrode or interface cable, hold the connector and pull it off. If you pull at the cable, it may cause a breakage.
- The touch panel is capacitance-type. Make sure to turn OFF the power before cleaning the panel. If you wipe it with the power ON, it may cause instrument malfunction.
- RS-232C or USB communication between the instrument and a personal computer may fail because of environmental conditions, such as (radio/electromagnetic) noise.

Environmental conditions for use and storage

- Temperature: 0°C to 45°C
- Humidity: under 80% in relative humidity and free from condensation

Avoid the following conditions:

- Dusty environment
- Strong vibration
- Direct sunlight
- Corrosive gas environment
- Close to an air-conditioner
- Direct wind

Transportation

When transporting the instrument, repackage it in the original package box. Otherwise, it may cause instrument breakage.

Disposal

Standard solution used for the calibration must be under neutralized before the disposal. As for the disposal of the meter, treat it as an industrial waste.

Description in This Manual

NOTE

This interprets the necessary points for correct operation and notifies the important points for handling the instrument.

REF

This indicates the part of where to refer the informatione.

— HINT!——

This indicates reference information.

CONTENTS

Chapter 1 About the F-74BW 1				
Chapter 2 Information of the F-74BW . 2				
2.1 Measurement Items 2				
2.2 Packing List 3				
2.3 Names and Functions 4 2.3.1 Names of Each Part 4 2.3.2 Display 5 2.3.3 Operation Key 6				
2.4 Measurement (MEAS) Flow Chart				
2.4.1 CH1				
2.5 Setting (SET) Flow Chart				
2.6 Data (DATA) Flow Chart				
Chapter 3 Basic Operation				
3.1 Preparation 11				
3.1.1 Assembling the Electrode Stand				
3.1.2 Electrode Connection				
3.1.3 Connecting the Power Source				
3.1.5 Connecting the Personal Computer				
3.2 Settings 16				
3.2.1 pH Standard Solution Setting163.2.2 Temperature Compensation Function Setting183.2.3 Auto Data Memory Setting203.2.4 pH Calibration Interval Setting21				
3.2.5 Sample ID Setting 22 3.2.6 ION Unit Setting 23				
3.2.7 ION Valency Setting				
3.2.8 Conductivity Unit Setting				

CONTENTS

	3.2.9 Conductivity Temperature Coefficient Setting	26
	3.2.10 Conductivity Cell Constant Setting	
	3.2.11 Salinity Unit Setting	
	3.2.12 Auto Printer Setting	
	3.2.13 Maintenance	
3.3	Calibration	
	3.3.1 pH Calibration	45
	3.3.2 ION Calibration	
	3.3.3 Conductivity Calibration	66
	3.3.4 Salinity Calibration	
	3.3.5 Calibration Value Printing	
3.4	Measurement 7	
	3.4.1 Measurement Value Hold	71
	3.4.2 Measurement Items Changing	72
	3.4.3 Auto Data Memory	
	3.4.4 Measurement Value Printing	
3.5	Data Operation 7	78
	3.5.1 Data Memory of Measurement Value	
	3.5.2 Memory Data Display	
	3.5.3 Memory Data Printing	
	3.5.4 Calibration History Display	
	3.5.5 Calibration History Printing	
	3.5.6 pH Check History Display	
	3.5.7 pH Check Data Printing	

4.3	CON		33
4.4	Τrοι	ubleshooting	100
	4.4.1	Error Message Chart	100
	4.4.2	More Troubleshooting	107

CONTENTS

Chapte	r 5 Specifications 11	13
5.1	Specifications 1	13
5.2	Default Settings 1	16
5.3	Options 1	17

Chapter 1 About the F-74BW

The desktop type pH meter F-74BW is optimized for laboratory measurement, and allows you to measure pH, ORP, ION, conductivity, resistivility, salinity, and TDS.

It provides a comfortable measurement environment with the design easy to wipe clean without roughness on the surface, and the free-standing electrode stand and the custom LCD for ease of measurement.

Chapter 2 Information of the F-74BW

2.1 Measurement Items

Measurement item	F-74BW
рН	
ORP (mV)	
Ion	
Conductivity	
Salinity	
Resistivity	
TDS (Total Dissolved Solids)	
Temperature	

2.2 Packing List


The F-74BW's package includes the following.


Electrodes are required for the measurement in addition to this package. If you want to buy an electrode, contact HORIBA sales representative.

Product	Quantity
Main unit	1
Electrode stand	1
Instruction Manual	1
Quick Manual	1
AC adapter	1
Ferrite core	1

2.3 Names and Functions


2.3.1 Names of Each Part





Connector for USB communication with a personal computer

2.3.2 Display



2.3.3 Operation Key

This instrument employs capacitance switches.

You can not operate them with thick gloves. Operate them with bare hands or thin rubber gloves.


Operation key		Function
CAL	CAL key	Starts/stops calibration or enters the calibration mode.
Ы		
MEAS	MEAS key	Starts/stops measurement or enters the measurement mode.
\sim		Applies the setting values of the setting mode.
DATA	DATA key	Enters the data operation mode.
F		Printouts calibration value.
SET	SET key	Enters the setting mode.
2		Cancels the setting values of the setting mode.
MODE	MODE key	Toggles the measurement items.
$\langle \rangle$		
\bigcirc	▲ key	Increases the values.
	▼ key	Decreases the values. Switches the setting item.
		Switches the CH items.
•	∢ ▶ key	Change number of digits, selects functions.
		Printouts measurement values. (key)
	ENTER key	Settlement, execution, prinout
		(except measurement values) Start/stops automatic data memory.
Ċ	POWER key	Turns ON or OFF the power.

Press the POWER key for 1 second or longer to turn ON the power, for 2 seconds or longer to turn OFF the power.


NOTE

2.4 Measurement (MEAS) Flow Chart

2.4.1 CH1




2.4.2 CH2




2.5 Setting (SET) Flow Chart

2.5.1 CH1






2.5.2 CH2



2.6 Data (DATA) Flow Chart




Chapter 3 Basic Operation


3.1 Preparation

3.1.1 Assembling the Electrode Stand

- **1.** Attach the stand shaft to the stand base.
- **2.** Attach the stopper and the stand arm to the stand shaft.



3.1.2 Electrode Connection




Electrode connector


1. Insert the groove of electrode plug by fitting it with the socket pin of the instrument.




If the pin and groove are misaligned, do not insert the plug with force.



2. Turn the electrode connector to the right along the groove to plug the connector.




3. Put the connector cover on the connector.



• Temperature connector


1. Insert the temperature connector into the jack socket on the instrument.



If the temperature connector is unconnected or the connection is wrong, the MTC set temperature is displayed as the sample temperature.

3.1.3 Connecting the Power Source

1. Insert the AC adapter cable by fitting with the connector socket of in the instrument.




- NOTE
- Do not insert the cable with force when the connectors do not match.
- Attach the provided ferrite core to the AC adapter cable.

3.1.4 Connecting and Setting the Printer

• Connecting the Printer

1. Insert the printer cable by fitting with the connector socket of the instrument.



The following printer is possible.

Printer

- CITIZEN CBM-910-24RJ120 V: plain paper type (Parts No.: 3014030146)
- CITIZEN CBM-910-24RJ230 V:
- plain paper type (Parts No.: 3014030147)
- Optional printer cable (Parts No.: 3014030148) is required.

NOTE


- · Make sure to use an appropriate cable for the printer.
- Make sure to power OFF the instrument before connecting a printer.
- When you do not connect a printer with the instrument, disconnect the printer cable and put the rubber cap firmly on the connector sockets on the instrument.

• Setting the Printer

Refer to the instruction manual of the printer for settings and operations of the printer.

- Set the DIP switch No. 6 to ON and No. 7 to OFF, and then set printer paper and ink ribbon. Keep the LF key held down.
- (2) Keep the SEL key held down. The printer prints output when the SEL key is being pressed.

3.1.5 Connecting the Personal Computer



Personal computer communication connector

•Use proprietary cables to connect with a personal computer.

Proprietary cable Parts name: USB cable (1 m) Parts No.: 3200373941

- Make sure that the transfer formats of the measuring instrument and personal computer are same. Otherwise, communication may fail due to a communication error or the online mode start failure. If you change the transfer formats, power OFF both of the instrument and the personal computer once, and then turn on them again.
- •For the details of communication commands, register with our website and see the free download page of manuals.

3.2 Settings

This section describes the procedures of the instrument condition settings for measurement, calibration, and maintenance.


The settings apply when you press the MEAS key to return to the measurement screen.

If you press the SET key during setting or after pressing the ENTER key, the settings are cancelled.

3.2.1 pH Standard Solution Setting


There are 3-type settings of standard solutions used for calibration; NIST, USA, and CUST. NIST is selected initially for pH standard solutions. This setting is possible only for CH1.

1. Press the SET key on the measurement screen is displayed to enter the pH standard solution setting screen.



2. Press the ENTER key to enter the pH standard solution change screen.

3. Press the ◀ or ► key to change NIST, USA, and CUST, and press the ENTER key.



specification solution (CUST) Bottle mark 1 2 3 4 5

NOTE

The calibration value of the pH7 standard solution differs between NIST and USA.

NIST specification: pH6.865 (at 25°C)

US specification: pH7.000 (at 25°C)


When you use NIST standard solution for calibration, make sure to set NIST here.

3.2.2 Temperature Compensation Function Setting


You can perform temperature compensation for pH electrode outputs.

This setting can be set to each of CH1 and CH2.

 Press the SET key on the measurement screen, and press the ▼ key to enter the temperature compensation function setting screen.



- **2.** Press the ENTER key to enter the temperature compensation function change screen.
- 3. Press the ◀ or ▶ key to change ATC or MTC, and press the ENTER key.



If you select MTC, press the < or
key to select a digit and press the ▲ or
key to change the value, and then press the ENTER key. Setting range: 0°C to 100°C



ATC

In Automatic Temperature Compensation (ATC), the instrument detects the solution temperature. When the temperature sensor is connected, the current temperature of the solutions is displayed automatically. If the sensor is not connected, the set value for MTC (default: 25°C) is displayed.

MTC


In Manual Temperature Compensation (MTC), the temperature sensor of the electrode is not used, and the solution temperature is entered manually.

3.2.3 Auto Data Memory Setting

Data memory can be performed automatically at a regular interval.

This setting applies for both of CH1 and CH2.

 Press the SET key on the measurement screen, and press the ▲ or ▼ key to enter the auto data memory setting screen.



- **2.** Press the ENTER key to enter the auto data memory change screen.
- Press the ◀ or ► key to change ON or OFF, and press the ENTER key.



If you select ON, press the < or
 key to select a digit and press the ▲ or
 key to change the value, and then press the ENTER key.

Setting range: 2 seconds to 24 hours

3.2.4 pH Calibration Interval Setting

This setting allows you to ensure calibration at a regular interval. When the set inteval (days) elapses after the previous calibration, ERROR No.0008 (Calibration cycle error) is displayed.

This setting is possible only for CH1.

1. Press the SET key on the measurement screen to enter the pH calibration interval setting screen.

- **2.** Press the ENTER key to enter the calibration interval change screen.
- Press the ◀ or ► key to change ON or OFF, and press the ENTER key.



If you select ON, press the < or
 key to select a digit and press the ▲ or
 key to change the value, and then press the ENTER key.

Setting range: 1 to 999 days

3.2.5 Sample ID Setting

If a sample ID is set here, it is recorded together with the measurement data at the data memory and it allows you to search data easily.

This setting applies for both of CH1 and CH2.

 Press the SET key on the measurement screen, and press the ▲ or ▼ key to enter the sample ID setting screen.



- **2.** Press the ENTER key to enter the sample ID change screen.
- Press the ◀ or ▶ key to select digit and press the ▲ or ▼ key to change value, and press the ENTER key.

Setting range: 0 to 9999




3.2.6 ION Unit Setting

Select the unit of ion measurement between g/L and mol/L.


If you use different unit settings between calibration and measurement, you can not obtain correct measurement results. When you change the unit setting, perform calibration again.

This setting is possible only for CH1.

 Press the SET key on the measurement screen, and press the ▲ or ▼ key to enter the ion unit setting screen.



- **2.** Press the ENTER key to enter the ion unit change screen.
- **3.** Select the unit to g/L or mol/L, and press the ENTER key.



3.2.7 ION Valency Setting

Select the valency of the ion to be measured.

If you use different ion valency settings between calibration and measurement, you can not obtain correct measurement results. When you change the ion valency setting, perform calibration again.

This setting is possible only for CH1.

 Press the SET key on the measurement screen, and press the ▲ or ▼ key to enter the ion valency setting screen.



- **2.** Press the ENTER key to enter the ion valency change screen.
- **3.** Select the ion valency, and press the ENTER key.



3.2.8 Conductivity Unit Setting


Select the unit of electric conductivity measurement among S/cm, S/m, and mS/cm FIX. When mS/cm FIX is selected, the unit is fixed to mS/cm.

This setting is possible only for CH2.

 Press the SET key on the measurement screen, and press the ▲ or ▼ key to enter the conductivity unit setting screen.



- **2.** Press the ENTER key to enter the conductivity unit change screen.
- Press the ◀ or ► key to select S/cm, S/m or mS/ cm FIX, and press the ENTER key.




3.2.9 Conductivity Temperature Coefficient Setting


Set the coefficient used for temperature compensation in electric conductivity measurement (Input range: 0.00% to 10.00%). The default is 2.00%.

This setting is possible only for CH2.


 Press the SET key on the measurement screen, and press the ▲ or ▼ key to enter the conductivity unit setting screen.



- **2.** Press the ENTER key to enter the conductivity temperature coefficient change screen.
- Press the ◀ or ► key to change ON or OFF, and press the ENTER key .



4. Press the < or
key to select digit and press the
▲ or
key to change value, and ENTER key.
Setting range: 0.00% to 10.00%



3.2.10 Conductivity Cell Constant Setting

Set the cell constant of electrode (Setting range: 0.700 to 1.300). When using a electrode for the instrument for the first time, set the cell constant written on the label of the electrode.

This setting is possible only for CH2.

 Press the SET key on the measurement screen, and press the ▲ or ▼ key to enter the conductivity cell constant setting screen.



2. Press the ENTER key to enter the conductivity cell constant change screen.

When the \blacktriangleleft or \blacktriangleright key is pressed, the display changes depending on the current unit setting.

For the SI unit (S/m)

 $\times 10 \rightarrow \times 100 \rightarrow \times 1000 \rightarrow \times 10$

For the conventional units (S/cm, mS/cm FIX)

 $\times 0.1 \rightarrow \times 1 \rightarrow \times 10$



3. Press the ENTER key.

4. Press the < or
 key to select digit and press the
 ▲ or
 ▲ key to change value, and press the ENTER key.

Setting range: 0.700 to 1.300




3.2.11 Salinity Unit Setting


You can select the unit of salinity measurement between PPT and %.

This setting is possible only for CH2.

 Press the SET key on the measurement screen, and press the ▲ or ▼ key to enter the salinity unit setting screen.



2. Press the ENTER key to enter the salinity unit change screen.



Press the ◀ or ► key to select PPT or %, and press the ENTER key.

HORIBA

3.2.12 Auto Printer Setting

You can set automatic printout when a printer is connected.


When this setting is ON, automatic printout is executed at the following occasions.

- You perform data memory during measurement.
- Calibration or check is completed in the calibration mode.

This setting applies for both of CH1 and CH2.

 Press the SET key on the measurement screen, and press the ▲ or ▼ key to enter the auto printer setting screen.

- **2.** Press the ENTER key to enter the auto printer setting change screen.
- Press the ◀ or ▶ key to select ON or OFF, and press the ENTER key.



Chapter 3 Basic Operation 3.2 Settings

3.2.13 Maintenance

The functions on the maintenance screen helps you to use the instrument in good conditions.

• LCD Check

You can check the LCD.

 Press the SET key on the measurement screen, and press the ▲ or ▼ key to enter the maintenance screen.




2. Press the ENTER key to enter the maintenance item selecting screen.



3. Press ENTER key to indicate the all LCD.



4. Press the ENTER key to return to the maintenance item selecting screen.




Printer Test

You can check the printer connection and printout operation.


 Press the SET key on the measurement screen, and press the ▲ or ▼ key to enter the maintenance screen.




2. Press the ENTER key to enter the maintenance item selecting screen.




3. Press the ▲ or ▼ key to enter the printer test screen.



4. Press the ENTER key to start the printing.



5. Press the ENTER key to return to the printer test screen.



Auto Power Off Setting

You can set ON/OFF of the auto power off function. When this function is ON and no key operations are performed for the set period, the instrument is powered OFF automatically.

 Press the SET key on the measurement screen, and press the ▲ or ▼ key to enter the maintenance screen.




2. Press the ENTER key to enter the maintenance item selecting screen.

Press the ▲ or ▼ key to enter the auto power off setting screen.



- **4.** Press the ENTER key to enter the auto power off change screen.
- Press the ◀ or ► key to change ON or OFF, and press the ENTER key.



6. Press the ◀ or ▶ key to select a digit and press the ▲ or ▼ key to change the value, and then press the ENTER key and press the ENTER key to return to the auto power off change screen.

Setting range: 1 to 30 minutes



Date Setting

Perform this operation when you use the instrument for the first time or replace the lithium batteries.

 Press the SET key on the measurement screen, and press the ▲ or ▼ key to enter the maintenance screen.



2. Press the ENTER key to enter the maintenance item selecting screen.



Press the ▲ or ▼ key to enter the date setting screen.



Press the ENTER key to enter the year setting screen. Press the ◄ or ► key to select a digit and press the ▲ or ▼ key to change the value, and then press the ENTER key.



Press the < or
 key to select a digit and press the ▲ or
 key to change the value for month and day setting, and then press the ENTER key.



6. Press the ◀ or ▶ key to select a digit and press the ▲ or ▼ key to change the value for hour and minute setting, and press the ENTER key to return to the date setting screen.




Memory Data Clearing

You can delete the memory data of the instrument. This operation is irreversible. We recommend printing out or saving data into a PC ("3.1.5 Connecting the Personal Computer" P.15) before this operation.

 Press the SET key on the measurement screen, and press the ▲ or ▼ key to enter the maintenance screen.


2. Press the ENTER key to enter the maintenance item selecting screen.



3. Press the ▲ or ▼ key to enter the memory data clear screen.




4. Press the ENTER key to enter the memory data clear (Yes/No) screen.



Press the ◀ or ► key to select Yes, and press ENTER key to delete the memory data.

When the deletion is completed, "End" is displayed.



6. Press the ENTER key to return to the memory data clear screen.

Calibration Data Clearing

You can delete the calibration data of pH, ion, conductivity, and salinity to initialize them. Note that this operation is irreversible.

 Press the SET key on the measurement screen, and press the ▲ or ▼ key to enter the maintenance screen.




2. Press the ENTER key to enter the maintenance item selecting screen.



Press the ▲ or ▼ key to enter the calibration data clear screen.




4. Press the ENTER key and the ◀ or ► key to select the measurement item to be deleted.



Press the ENTER key to enter the calibration data deleting (Yes/No) screen.

(*Case of pH)



6. Press the ◀ or ▶ key to select Yes, and press ENTER key to delete the calibration data.
 When the deletion is completed, "End" is displayed.

7. Press the ENTER key to return to the calibration data clear screen.


Number of Memory Data Display

The number of memory data saved in the instrument is displayed. You can store up to 999 data.


 Press the SET key on the measurement screen, and press the ▲ or ▼ key to enter the maintenance screen.



2. Press the ENTER key to enter the maintenance item selecting screen.




Press the ▲ or ▼ key to enter the number of memory data display screen.



4. Press the ENTER key to display the number of the memory data currently saved in the instrument.



5. Press the ENTER key to return to the number of memory data display screen.



Temperature Calibration

You can perform calibration of the temperature sensor of an electrode.


This setting is possible for both CH1 and the CH2.

 Press the SET key on the measurement screen, and press the ▲ or ▼ key to enter the maintenance screen.



2. Press the ENTER key to enter the maintenance item selecting screen.

3. Press the ▲ or ▼ key to enter the temperature calibration screen.



- **4.** Press the ENTER key to enter the temperature display value change screen.
- Press the ◀ or ► key to select a digit and press the ▲ or ▼ key to change the value, and then press the ENTER key to return to the temperature calibration screen.



Setting Initialization


The settings are initialized.

Refer to "5.2 Default Settings" P.116 for default settings.


 Press the SET key on the measurement screen, and press the ▲ or ▼ key to enter the maintenance screen.



2. Press the ENTER key to enter the maintenance item selecting screen.



Press the ▲ or ▼ key to enter the initializing settings screen.




4. Press the ENTER key to enter the initializing settings (Yes/No) screen.



Press the ◀ or ▶ key to select Yes, and press ENTER key to initialize the settings.

The initializing is complete when "End" is displayed.




6. Press the ENTER key to return to the initializing settings screen.

Password Setting


This function prevents an unauthorized person from tampering the date of measurement data, or operating accidentally to change settings, delete memory data, etc.

 Press the SET key on the measurement screen, and press the ▲ or ▼ key to enter the maintenance screen.


2. Press the ENTER key to enter the maintenance item selecting screen.



Press the ▲ or ▼ key to enter the password settings screen.



- **4.** Press the ENTER key to enter the password change screen.
- Press the ◀ or ► key to change ON or OFF, and press the ENTER key.



6. If you select ON, press the ◀ or ▶ key to select a digit and press the ▲ or ▼ key to change the value, and then press the ENTER key and press the ENTER key to return to the password change screen.

Setting range: 0000 to 9999

NOTE



If this setting is ON, password entry is required when you power ON the instrument. Enter the password and press the ENTER key.

If you forget the password, contact HORIBA sales representative.

3.3.1 pH Calibration

Up to 5 points of pH calibration is possible.

Repeatability Check of Calibration

After calibration is completed, immerse the electrode in the standard solution of pH7 and press the ENTER key to perform the repeatability check of calibration. It is necessary to calibrate the standard solution of pH7 in NIST or USA of the standard solution setting. This operation is valid only once after calibration.

Calibration with NIST Standard Solution

It is possible to calibrate pH2, pH4, pH7, pH9, and pH12 using NIST standard solutions.

In this manual, an example of 2 points calibration with the standard solutions of pH7 and pH4, which are used generally, is explained.

The calibration order is arbitrary. In this example, the procedures of calibration of pH7 are explained, and followed by ones of pH4, but you can perform calibration in an order as you like.

- **1.** Prepare the pH electrode. Wash the pH electrode with pure water (deionized), and wipe it off by filter paper or tissue paper.
- **2.** Open the filler port of the internal solution. During the calibration, the filler port should be kept open.
- **3.** Immerse the pH electrode in the standard solution of pH7 to calibrate the first point.

4. Press the ▲ or ▼ key on the measurement screen to light up ACT for CH1.



5. Press the CAL key to enter the NIST calibration screen.



6. Wait for the instantaneous value to be stable and press the CAL key.

The HOLD indicator blinks until the reading stabilizes. When the reading stabilizes, the HOLD indicator lights up and the reading is held.



The **7** mark, which means that the pH7 calibration is completed, is displayed on the screen.

• To stop calibration.

Press the CAL key while the HOLD indicator blinks.

• To redo calibration.

Press the CAL key again after the HOLD indicator lights up.

- 7. Wash the pH electrode with pure water (deionized), and wipe it off by filter paper or tissue paper.
- **8.** Immerse the pH electrode in the standard solution of pH4 to calibrate the second point.
- **9.** Press the CAL key to display the instantaneous value. Wait for the instantaneous value to be stable and press the CAL key.

The HOLD indicator blinks until the reading stabilizes. When the reading stabilizes, the HOLD indicator lights up and the reading is held.



The **4** mark, which means that the pH4 calibration is completed, is displayed on the screen.

10. To perform repeatability check of the calibration, immerse the pH electrode in the standard solution of pH7 and press the Enter key twice.

The measurement value of the pH7 standard solution is displayed on the screen.

11. Press the MEAS key to return to the pH measurement mode.

NOTE

You can redo the calibration procedures of each standard solution while you in the same calibration mode. But if you go back to the measurement mode once and return to the calibration mode again, redoing calibration will erase the previous set of calibration. You have to perform the entire calibration procedures from the beginning in this case.

Electrode status

You can check the electrode status after calibration.

Display	Description		
,ERR No display	Good condition. Electrode sensitivity: 93% to 100%		
Blinks	 Electrode sensitivity: 90% to 93% Check that the standard solutions are correct. Wash the electrode. " Washing electrodes" P.91 		
Lights up ERROR No.0004	 Electrode sensitivity: 85% to 90% " ● ERROR No.0005 Span calibration ERROR (pH)" P.103 It is the asymmetry potential error. " ● ERROR No.0004 Zero calibration 		
ERROR No.0005	ERROR" P.102 It is the sensitivity error. " ● ERROR No.0005 Span calibration ERROR (pH)" P.103		

Calibration with USA Standard Solution

It is possible to calibrate pH2, pH4, pH7, pH10, and pH12 using USA standard solutions.

In this manual, an example of 2 points calibration with the standard solutions of pH7 and pH4, which are used generally, is explained.


The calibration order is arbitrary. In this example, the procedures of calibration of pH7 are explained, and followed by ones of pH4, but you can perform calibration in an order as you like.

1. Prepare the pH electrode. Wash the pH electrode with pure water (deionized), and wipe it off by filter paper or tissue paper.

- **2.** Open the filler port of the internal solution. During the calibration, the filler port should be kept open.
- **3.** Immerse the pH electrode in the standard solution of pH7 to calibrate the first point.
- **4.** Press the ▲ or ▼ key on the measurement screen is displayed to light up ACT for CH1.



5. Press the CAL key to enter the USA calibration screen.



6. Check the stable indication value, and press the CAL key.

The HOLD indicator blinks until the reading stabilizes. When the reading stabilizes, the HOLD indicator lights up and the reading is held.



The **7** mark, which means that the pH7 calibration is completed, is displayed on the screen.

 To stop calibration: Press the CAL key while the HOLD indicator blinks.

- To redo calibration: Press the CAL key again after the HOLD indicator lights up.
- 7. Wash the pH electrode with pure water (deionized), and wipe it off by filter paper or tissue paper.
- **8.** Immerse the pH electrode in the standard solution of pH4 to calibrate the second point.
- **9.** Press the CAL key to display the instantaneous value. Wait for the instantaneous value to be stable and press the CAL key.

The HOLD indicator blinks until the reading stabilizes. When the reading stabilizes, the HOLD indicator lights up and the reading is held.



The **4** mark, which means that the pH4 calibration is completed, is displayed on the screen.

10. To perform repeatability check of the calibration, immerse the pH electrode in the standard solution of pH7 and press the Enter key twice.

The measurement value of the pH7 standard solution is displayed on the screen.

11. Press the MEAS key to return to the pH measurement mode.

NOTE

You can redo the calibration procedures of each standard solution while the calibration mode is continued. But if you go back to the measurement mode once and return to the calibration mode again, the calibration data are all cleared. You have to perform the entire calibration procedures from the beginning in this case.

Electrode status

You can check the electrode status after calibration.

Display	Description		
,ERR No display	Good condition. Electrode sensitivity: 93% to 100%		
Blinks	 Electrode sensitivity: 90% to 93% Check that the standard solutions are correct. Wash the electrode. " ● Washing electrodes" P.91 		
Lights up	Electrode sensitivity: 85% to 90% " ● ERROR No.0005 Span calibration ERROR (pH)" P.103		
ERROR No.0004	It is the asymmetry potential error. " ● ERROR No.0004 Zero calibration ERROR" P.102		
ERROR No.0005	It is the sensitivity error. " ● ERROR No.0005 Span calibration ERROR (pH)" P.103		

Calibration with CUST Screen


Up to 5 points calibration can be performed using desired standard solutions.

- **1.** Prepare the pH electrode. Wash the pH electrode with pure water (deionized), and wipe it off by filter paper or tissue paper.
- **2.** Open the filler port of the internal solution. During the calibration, the filler port should be kept open.
- **3.** Immerse the pH electrode in the desired standard solution to calibrate the first point.

4. Press the ▲ or ▼ key on the measurement screen to light up ACT for CH1.



5. Press the CAL key to enter the CUST calibration screen.



- 6. Press the ◄ or ▶ key to select digit and press the
 ▲ or ▼ key to change value, to set the value of the desired standard solution.
- **7.** Check the stable indication value, and press the CAL key.

The HOLD indicator blinks until the reading stabilizes. When the reading stabilizes, the HOLD indicator lights up and the reading is held.




The **mark**, which means that the first point calibration is completed, is displayed on the screen.

- To stop calibration: Press the CAL key while the HOLD indicator blinks.
- To redo calibration:

Press the CAL key again after the HOLD indicator lights up.

- **8.** Wash the pH electrode with pure water (deionized), and wipe it off by filter paper or tissue paper.
- **9.** Immerse the pH electrode in the standard solution of pH4 to calibrate the second point.
- **10.**Press the ▲ or ▼ key to set the value of the random standard solution.
- **11.** Press the CAL key to display the instantaneous value. Wait for the instantaneous value to be stable and press the CAL key.

The HOLD indicator blinks until the reading stabilizes. When the reading stabilizes, the HOLD indicator lights up and the reading is held.



The **2** mark, which means that the second point calibration is completed, is displayed on the screen.

- **12.**As for calibration of the third and later points, perform the calibration from the step 7.
- **13.**When the calibration is completed, press the MEAS key to return to the measurement mode.

3.3.2 ION Calibration

Ion electrodes can measure ion concentrations. But, a skill is necessary for more accurate measurement. Refer to the instruction manual of the electrode and the following descriptions.

Up to 5 points ion calibration is possible.

1. Addition of ionic strength conditioner to sample

The effect of the ionic strength may cause an error in measurement values. Adjust ionic strength by adding the designated reagent to the sample and standard solutions. And it will allow you to make accurate measurement.

2. pH effect of sample

The applicable pH range depends on the type and structure of the ion electrode. Check that the pH of each sample is within the applicable range. If the pH is out of the range, adjust the sample pH using the reagent that does not contain any ions to be measured or interfering.

3. Electric potential slope measured by ion measurement method

Ion concentration measurement is affected by the temperature of solution to be measured in accordance with the Nerst formula. Also, no stirring will cause slow response, measurement failure in low concentrations, or the fluctuation of readings. Use a temperature controlled bath with a stirrer for standard solution calibration and sample measurement.

4. Measurement of ion concentration

When a particular ion exist within the solution to be measured, the responsive membrane of an ion electrode generates an electric potential corresponding to the concentrations of the ion. The generated potential can be measured using an ion meter as the potential difference from the reference electrode.

With ion electrodes, the measured potential is generally proportional to logarithm of the ion activity in the measured solution, and expressed in the following formula.


 $E=E_0+(2.303RT/nF)\log \gamma C$

- E: Measured electric potential (V)
- E₀: Standard potential (V), determined according to the system. This includes the standard potential of the reference electrode and the liquid junction potential.
- F: Faraday constant (96,480 Cmol⁻¹)
- R: General gas constant (8.314 JK⁻¹mol⁻¹)
- T: Absolute temperature (K)
- n: Ion charge
- γ : Activity coefficient
- C: Ion concentration (mol/L)

The above formula is called "Nernst's equation" and is the basis for measuring ion concentrations using an ion electrode.

(2.303RT/nF), which is called the potential slope, incline, slope, or Nernst's factor, is the change in potential generated when the ion concentrations change by a factor of 10. This change in potential is called the potential slope, incline, slope, or Nernst's factor. According to the equation, if you determine the values of potential slope and E0 by calibrating the standard solutions, you can obtain the concentrations of the target ion by obtaining the potential difference, E of the ion electrode inside the measured solution.

In actual measurement, an ion electrode measures the ion concentrations, therefore, a linear relationship forms between the logarithm of ion concentrations and the electrode potential. Quantitative analysis using an ion electrode requires either an ion meter that has an antilog calculation function, or the creation of a calibration curve using semilog graph paper.



Calibration curve for univalent positive ion electrode



5. Interfering ions effect of sample

If sample measurement is not successful even after the countermeasures of 1 through 3 mentioned above, the sample may contain interfering ions.

How to reduce the effect of the interfering ion, contact HORIBA sales representative.

When (coexistence limit value) \times (measurement ion concentration) is larger than (coexistence ion concentration), the effect of the interfering ion becomes small.

Harmful chemicals

Some ion electrodes are used with hazardous standard solutions. Handle them with care.

If the internal solution comes in contact with the skin, wash it off immediately. If it gets into eyes, flush with plenty of water and then consult a doctor.

Information on Applicable Electrodes

(1) Ion electrode (single-pole type)

Electrode type	Model and ion types	 Measurement range pH allowable range*¹ Electric potential slope*² Temperature range 	Interfering ion Coexistence allowable limit ^{*3}	Applicable reference electrode and its internal liquid type ^{*4}
Solid mem- brane	8001- 10C Cyanide (CN⁻)	 0.03 mg/L to 2600 mg/L CN⁻ (10⁻¹ to 10⁻⁶ mol/L CN⁻) pH12 to pH13 (2.6 mg/L CN⁻¹) -59 mV 0°C to 50°C 	Cl ⁻ , F ⁻ , Br ⁻ , NO ₃ ⁻ , HCO ₃ ⁻ , SO ₄ ²⁻ , SO ₃ ²⁻ =10000, S ₂ O ₃ ²⁻ =1 l ⁻ =0.1 S ²⁻ , MnO ₄ ⁻ = Not applicable (At 10 ⁻³ mol/L CN ⁻)	2060A-10T, 2565A-10T 3.33 mol/L KCI (250 g/L potassium chloride)
	8002- 10C Chloride (CI⁻)	 0.4 mg/L to 35000 mg/L Cl⁻ (1 to 10⁻⁵ mol/L Cl⁻) pH3 to pH11 (350 mg/L Cl⁻) -59 mV 0°C to 50°C 	NO ₃ ⁻ , F ⁻ , HCO ₃ ⁻ , SO ₄ ²⁻ , PO ₄ ²⁻ =1000 SCN ⁻ =0.3 MnO ₄ ⁻ =0.1 Br ⁻ =0.03 S ₂ O ₃ ²⁻ , S ²⁻ , I ⁻ , Ag ⁺ , Hg ²⁺ = Not applicable (At 10 ⁻³ mol/L Cl ⁻)	2565A-10T 1 mol/L KNO ₃ (100 g/L Potassium nitrate)
	8003- 10C Sulfide (S ²⁻)	 0.3 mg/L to 32000 mg/L S²⁻ (1 to 10⁻⁵ mol/L S²⁻) pH12 to pH14 (3.2 mg/L S²⁻) -29 mV 0°C to 50°C 	NO ₃ ⁻ , Br ⁻ =10000 I ⁻ , F ⁻ , CI ⁻ , PO ₄ ²⁻ , SO ₄ ²⁻ =1000 S ₂ O ₃ ²⁻ =10 CN ⁻ = Not applicable (At 10 ⁻⁴ mol/L S ⁻)	2060A-10T, 2565A-10T 3.33 mol/L KCI (250 g/L potassium chloride)
	8004- 10C Iodide (I ⁻)	 0.01 mg/L to 13000 mg/L l⁻ (10⁻¹ to 10⁻⁷ mol/L l⁻) pH2 to pH11 (1300 mg/L l⁻) -59 mV 0°C to 50°C 	Cl ⁻ , F ⁻ , PO ₄ ²⁻ , SO ₄ ²⁻ =10000, Br ⁻ =1000 NO ₃ ⁻ =100 $S_2O_3^{2-}=10$ MnO ₄ ^{-,} S ^{2-,} CN ⁻ = Not applicable (At 10 ⁻³ mol/L l ⁻)	2060A-10T, 2565A-10T 3.33 mol/L KCI (250 g/L potassium chloride)

Chapter 3 Basic Operation 3.3 Calibration

Electrode type	Model and ion types	 Measurement range pH allowable range*¹ Electric potential slope*² Temperature range 	Interfering ion Coexistence allowable limit ^{*3}	Applicable reference electrode and its internal liquid type ^{*4}
Solid mem- brane	8005- 10C Bromide (Br ⁻)	 0.8 mg/L to 80000 mg/L Br⁻ (1 to 10⁻⁵ mol/L Br⁻) pH1.5 to pH11.5 (800 mg/L Br⁻) -59 mV 0°C to 50°C 	F ⁻ , NO ₃ ⁻ , SO ₄ ²⁻ =1000 Cl ⁻ , PO ₄ ²⁻ =100 MnO ₄ ⁻ =1 S ₂ O ₃ ²⁻ , S ²⁻ , l ⁻ , CN ⁻ = Not applicable (At 10 ⁻³ mol/L Br ⁻)	2065A-10T 1 mol/L KNO ₃ (100 g/L Potassium nitrate)
	8006- 10C Copper (Cu ²⁺)	 0.06 mg/L to 6400 mg/L Cu²⁺ (10⁻¹ to 10⁻⁶ mol/L Cu²⁺) pH2 to pH6 (6.4 mg/L Cu²⁺) 29 mV 0°C to 50°C 	Zn ²⁺ , Mg ²⁺ , NH ₄ ⁺ =10000 Ni ²⁺ , Na ⁺ =1000 Fe ²⁺ =0.1	2060A-10T, 2565A-10T 3.33 mol/L KCI (250 g/L potassium chloride)
	8007- 10C Cadmium (Cd ²⁺)	 0.1 mg/L to 11000 mg/L Cd²⁺ (10⁻¹ to 10⁻⁶ mol/L Cd²⁺) pH3 to pH8 (11 mg/L Cd²⁺) 29 mV 0°C to 50°C 	K ⁺ , Mg ²⁺ =100000 Zn ²⁺ , NH ₄ ⁺ =10000 Ni ²⁺ =1000 Cr ³⁺ , Fe ²⁺ =100 Fe ³⁺ =1 Pb ²⁺ =0.1 Cu ²⁺ , Hg ²⁺ , Ag ²⁺ = Not applicable	2060A-10T, 2565A-10T 3.33 mol/L KCI (250 g/L potassium chloride)
	8008- 10C Lead (Pb ²⁺)	 2.0 mg/L to 20000 mg/L Pb²⁺ (10⁻¹ to 10⁻⁵ mol/L Pb²⁺) pH4.5 to pH6.5 (20 mg/L Pb²⁺) 29 mV 0°C to 50°C 	NH_4^+ , $K^+=1000$ Ni^{2+} , Mg^{2+} , $Zn^{2+}=100$ $Cd^{2+}=10$ $Cr^{3+}=1$ $Fe^{3+}=0.01$ Cu^{2+} , Hg^{2+} , Ag^+ , $S^{2-}=$ Not applicable	2565A-10T 1 mol/L KNO ₃ (100 g/L Potassium nitrate)
	8009- 10C Thiocya nate (SCN ⁻)	 ① 0.6 mg/L to 5800 mg/ L SCN⁻ (10⁻¹ to 10⁻⁵ mol/L SCN⁻) ② pH2 to pH12 (5.8 mg/L SCN⁻) ③ -59 mV ④ 0°C to 50°C 	NO_3^- , CIO_4^- , $SO_4^{-2-}=10000$ $CI^-=100$ $Br^-=1$ $S_2O_3^{-2-}$, S^{-2-} , I^- , $CN^-=$ Not applicable (At 10^{-3} mol/L SCN ⁻)	2565A-10T 1 mol/L KNO ₃ (100 g/L Potassium nitrate)

Chapter 3 Basic Operation 3.3 Calibration

Electrode type Solid mem- brane	Model and ion types 8010- 10C Fluoride	 Measurement range pH allowable range*¹ Electric potential slope*² Temperature range 0.02 mg/L to 19000 mg/L F⁻ (1 to 10⁻⁶ mol/L F⁻) 	Interfering ion Coexistence allowable limit ^{*3} OH ⁻ =10 Disturbance may occur due to the formation of	Applicable reference electrode and its internal liquid type ^{*4} 2060A-10T, 2565A-10T 3.33 mol/L KCI
	(F ⁻)	 ② pH4 to pH10 (19 mg/L F⁻) ③ -59 mV ④ 0°C to 50°C 	complex when polyvalent ions such as Al ³⁺ , Fe ³⁺ coexist.	(250 g/L potassium chloride)
	8011- 10C Silver (Ag ⁺)	 0.01 mg/L to 110000 mg/L Ag⁺ (1 to 10⁻⁷ mol/L Ag⁺) pH2 to pH10 (1 mg/L Ag⁺) 59 mV 0°C to 50°C 	Cu ²⁺ , Cd ²⁺ , Pb ²⁺ , Zn ²⁺ , Mg ²⁺ , Ca ²⁺ , Na ⁺ , K+=1000 Hg ²⁺ = Not applicable	2565A-10T 1 mol/L KNO ₃ (100 g/L Potassium nitrate)
Dia- phragm	5002A- 10C Ammonia (NH ₃)	 0.1 mg/L to 1000 mg/ L NH₃ Adjust to pH12 or more. -59 mV 0°C to 50°C 	Substances forming acidic or basic gases in alkali solutions of pH12 or higher.(volatile amines, etc.)	Because of the combination type ion electrode, the reference electrode is not reqiured.
Glass	1512A- 10C Sodium (Na ⁺)	 2.3 mg/L to 23000 mg/L Na⁺ (10 to 10⁻⁴ mol/L Na⁺) more than pH4.5 (2300 mg/L Na⁺) +59 mV 0°C to 60°C 	Ca ²⁺ =500 NH ₄ ⁺ =20 K ⁺ , Li ⁺ =10	2565A-10T 0.1 mol/L KCI (7.5 g/L potassium chlo- ride) or Sample
Liquid mem- brane	8201- 10C Nitrate (NO ₃ ⁻)	 0.62 mg/L to 62000 mg/L NO₃⁻ (1 to 10⁻⁵ mol/L NO₃⁻) pH3 to pH7 (62 mg/L NO₃⁻) -55 mV 0°C to 50°C 	SO ₄ ²⁺ =1000 以上 CH ₃ COO ⁻ =300 F ⁻ =200 CI ⁻ =40 NO ₂ ⁻ =3 Br ⁻ =2 I ⁻ =0.1 CIO ₄ ⁻ =0.03 (At 10 ⁻³ mol/L NO ₃ ⁻)	2565A-10T 1 mol/L KCI (75 g/L potassium chloride)

Chapter 3 Basic Operation 3.3 Calibration

Electrode type	Model and ion types	 Measurement range pH allowable range*¹ Electric potential slope*² Temperature range 	Interfering ion Coexistence allowable limit ^{*3}	Applicable reference electrode and its internal liquid type ^{*4}
Liquid mem- brane	8202- 10C Potassium (K ⁺)	 0.04 mg/L to 39000 mg/L K⁺ (1 to 10⁻⁶ mol/L K⁺) pH5 to pH11 (3.9 mg/L K⁺) +58 mV 0°C to 50°C 	Li ⁺ , Na ⁺ , Mg ²⁺ , Ca ²⁺ , Sr ²⁺ , Ba ²⁺ = more than 1000 NH ₄ ⁺ =70 Cs ⁺ =3 Rb ⁺ =0.4 (At 10 ⁻⁴ mol/L K ⁺)	2565A-10T 3.33 mol/L NaCl (200 g/L Sodium chloride)
	8203- 10C Calcium (Ca ²⁺)	 0.40 mg/L to 40080 mg/L Ca²⁺ (1 to 10⁻⁵ mol/L Ca²⁺) pH5 to pH11 (4.0 mg/L Ca²⁺) +29 mV 0°C to 50°C 	Na ⁺ , K ⁺ , Ba ⁺ , NH ₄ ⁺ = more than 1000 Mg ²⁺ =1000 Mn ²⁺ =500 Co ²⁺ =350 Ni ²⁺ , Cu ²⁺ =70 Sr ²⁺ =50 Fe ²⁺ , Zn ²⁺ =1 Fe ³⁺ =0.1 (At 10 ⁻⁴ mol/L Ca ²⁺)	2060A-10T 2565A-10T 3.33 mol/L KCI (250 g/L potassium chloride)

(2) Combination type ION electrodes

The combination type ion electrode are designed to allow its chip type ion electrode to be replaced easily. If you measure some of the ion species listed in the following table in addition to those included in the standard specifications, change the chip type ion electrode and the internal liquid in the reference part of the combination type ion electrode. These changes allows you to perform measurement un up to 5 ion species.

Each chip type ion electrode is available for purchase it is provided as one of the standard accessories. In accordance with the target ion species, select appropriate chip electrodes from the following table. The applicable temperature range is 0°C to 50°C for all the electrodes.

Model and ion type	Applica- ble chip electrode type (quantity)	 Measurable range pH allowable range^{*1} Electric potential slope^{*2} 	Interfering ions and allowable limit for their coexistence ^{*3}	Reference electrode (outer cylinder) and internal liquid type ^{*4}
6560-10C Chloride Cl ⁻	7660 (1 pc)	 0.4 mg/L to 35000 mg/L Cl⁻ (1 to 10⁻⁵ mol/L Cl⁻) pH3 to pH11 (350 mg/L Cl⁻) -59 mV 	NO ₃ ⁻ , F ⁻ , HCO ₃ ⁻ , SO ₄ ²⁻ , PO ₄ ²⁻ = more than 1000 SCN ⁻ =0.3 MnO ₄ ⁻ =0.1 Br ⁻ =0.03 S ₂ O ₃ ²⁻ , S ²⁻ , I ⁻ , Ag ⁺ , Hg ²⁺ = Not applicable (At 10 ⁻³ mol/L Cl ⁻)	1 mol/L KNO ₃ (100 g/L Potassium nitrate)
6561-10C Fluoride F ⁻	7661 (1 pc)	 ① 0.02 mg/L to 19000 mg/ L F⁻ (1 to 10⁻⁶ mol/L F⁻) ② pH4 to pH10 (19 mg/L F⁻) ③ -59 mV 	OH ⁻ =10 Disturbance may occur due to the formation of complex when polyvalent ions such as Al ³⁺ , Fe ³⁺ coexist.	3.33 mol/L KCI (250 g/L potassium chloride)
6581-10C Nitrate NO ₃ ⁻	7681 (2 pcs)	 0.62 mg/L to 62000 mg/ L NO₃⁻ (1 to 10⁻⁵ mol/L NO₃⁻) pH3 to pH7 (62 mg/L NO₃⁻) -55 mV 	SO_4^{2+} = more than 1000 $CH_3COO^-=300$ $F^-=200$ $CI^-=40$ $NO_2^-=3$ $Br^-=2$ $I^-=0.1$ $CIO_4^-=0.03$ (At 10^{-3} mol/L NO_3^-)	1 mol/L KCI (100 g/L potassium chloride)

Model and ion type	Applica- ble chip electrode type (quantity)	 Measurable range pH allowable range^{*1} Electric potential slope^{*2} 	Interfering ions and allowable limit for their coexistence ^{*3}	Reference electrode (outer cylinder) and internal liquid type ^{*4}
6582-10C Potassium K ⁺	7682 (2 pcs)	 0.04 mg/L to 39000 mg/ L K⁺ (1 to 10⁻⁶ mol/L K⁺) pH5 to pH11 (3.9 mg/L K⁺) 58 mV 		3.33 mol/L NaCl (200 g/L Sodium chloride)
6583-10C Calcium Ca ²⁺	7683 (2 pcs)	 0.40 mg/L to 40080 mg/ L Ca²⁺ (1 to 10⁻⁵ mol/L Ca²⁺) pH5 to pH11 (4.0 mg/L Ca²⁺) 29 mV 	Na ⁺ , K ⁺ , Ba ⁺ , NH ₄ ⁺ = more than 1000 Mg ²⁺ =1000 Mn ²⁺ =500 Co ²⁺ =350 Ni ²⁺ , Cu ²⁺ =70 Sr ²⁺ =50 Fe ²⁺ , Zn ²⁺ =1 Fe ³⁺ =0.1 (At 10 ⁻⁴ mol/L Ca ²⁺)	3.33 mol/L KCI (250 g/L potassium chloride)

- *1: The allowable pH range for the value specified in the parentheses under the measurable ion concentration. (pH range slightly changes depending on the ion concentrations.)
- *2: The change in electric potential at the electrode which occurs when the ion concentration increased by 10 times (25°C).
- *3: When the concentration of interfering ions has exceeds the value which is obtained by multiplying the concentrations of the measured ions by the allowable coexistence limit, the electric potential at the electrode is affected by interfering ions and this makes measurement impossible.

(The effect by the coexistence of interfering ions decreases as the allowable coexistence limit becomes larger or the concentrations of the measured ion become higher.)

*4: The internal liquid types are not limited to the solutions listed in the table. You may use other internal liquid so far as you have determined that its bleeding rarely interferes or affects the electric potential at the electrode and less frequently causes electric potential differences between the solutions.

- Prepare the standard solution. Refer to "Electrode Instruction Manual" for the preparative method.
- **2.** Calibrate the first point. Wash the ion electrode with pure water (deionized), and wipe it off by filter paper or tissue paper.
- **3.** Immerse the ion electrode in the calibration standard solution of the first point.
- **4.** Press the ▲ or ▼ key on the measurement screen to light up ACT for CH1.



- **5.** Press the MODE key on the measurement screen to enter the ion measurement screen.
- **6.** Press the CAL key on the ion measurement screen to enter the calibration screen.
- 7. Press the ◄ or ► key to select digit and press the
 ▲ or ▼ key to change value, to set the value of the desired standard solution.
- **8.** Wait for the instantaneous value to be stable and press the CAL key.

The HOLD indicator blinks until the reading stabilizes. When the reading stabilizes, the HOLD indicator lights up and the reading is held.



The **mark**, which means that the first point calibration is completed, is displayed on the screen.

- To stop calibration: Press the CAL key while the HOLD indicator blinks.
- To redo calibration: Press the CAL key again after the HOLD indicator lights up.
- **9.** Wash the ion electrode with pure water (deionized), and wipe it off by filter paper or tissue paper.
- **10.** Immerse the ion electrode in the calibration standard solution of the second point.
- **11.** Press the CAL key on the ion measurement screen to enter the calibration screen.
- 12.Press the < or
 key to select digit and press the
 ▲ or
 key to change value, to set the value of the desired standard solution.
- **13.**Press the CAL key to display the instantaneous value. Wait for the instantaneous value to be stable and press the CAL key.

The HOLD indicator blinks until the reading stabilizes. When the reading stabilizes, the HOLD indicator lights up and the reading is held.



The **2** mark, which means that the second point calibration is completed, is displayed on the screen.

• To stop calibration:

Press the CAL key while the HOLD indicator blinks.

• To redo calibration:

Press the CAL key again after the HOLD indicator lights up.

14.As for calibration of the third to fifth point, perform the calibration from step 8.

3.3.3 Conductivity Calibration

The cell constant of the electrode may change according to the using conditions. The recommended calibration of the cell constant is approximately once a year. The calibration updates the cell constant to be appropriate for the current electrode.

- **1.** Prepare the potassium chloride standard solution. Refer to "JIS K 0400" for the preparative method.
- **2.** Wash the electrode with pure water (deionized), and wipe it off by filter paper or tissue paper.
- **3.** Immerse the electrode in the potassium chloride standard solution of the second point.
- **4.** Press the ▲ or ▼ key when the measurement screen is displayed to light up ACT for CH2.



- **5.** Press the CAL key on the measurement screen is displayed to enter the calibration screen.
- 6. Press the ◄ or ► key to select digit and press the
 ▲ or ▼ key to change value, to set the value of the desired standard solution.

7. Check the stable indication value, and press the CAL key.

The HOLD indicator blinks until the reading stabilizes. When the reading stabilizes, the HOLD indicator lights up and the reading is held.



8. After the calibration is completed, press the MEAS key to return to the measurement screen.

3.3.4 Salinity Calibration

You can perform 1 point calibration using a salinity standard solution.

- Prepare the electrode. Wash the electrode with pure water (deionized), and wipe it off by filter paper or tissue paper.
- **2.** Immerse the electrode in the desired standard solution used for calibration.
- 3. Press the ▲ or ▼ key when the measurement screen is displayed to light up ACT for CH2.



- **4.** Press the MODE key on the measurement screen to enter the salt measurement screen.
- **5.** Press the CAL key to enter the salinity calibration screen.

- 6. Press the ◄ or ▶ key to select digit and press the
 ▲ or ▼ key to change value, to set the value of the desired standard solution.
- 7. Check the stable indication value, and press the CAL key.

The HOLD indicator blinks until the reading stabilizes. When the reading stabilizes, the HOLD indicator lights up and the reading is held.




8. After the calibration is completed, press the MEAS key to return to the measurement screen.

3.3.5 Calibration Value Printing

You can print out the data during calibration HOLD or calibration error.


- **1.** Check the printer settings and that the printer is connected correctly. (Refer to "3.1.4 Connecting and Setting the Printer" P.14.)
- **2.** Press the DATA key during calibration HOLD or calibration error. The currently displayed data is printed out. Printing examples are shown below.

• Printing at pH calibration



- ① Instrument model
- (2) Serial number
- 3 Electrode status
- (4) Asymmetry potential
- (5) Sensitivity
- 6 Calibration data

Printing at ION calibration




- 1 Instrument model
- 2 Serial number
- 3 Electrode status
- (4) Asymmetry potential
- 5 Sensitivity
- 6 Calibration data

Printing at COND calibration

- Inst. model :F-74BW-(1) <u>(</u>2) Inst. SN :123456789 CELL 3 0.973x1000 m-1-Calibration data :2011/03/10 Date Time :18:44 (4) 98.64 S/m 25.0°C MTC Temp. Coef : 2.00%/ C
- ① Instrument model
- (2) Serial number
- ③ Cell constant
- (4) Calibration data

Printing at salinity calibration



- ① Instrument model
- (2) Serial number
- 3 Calibration coefficient
- (4) Calibration data

3.4.1 Measurement Value Hold

The instrument can hold the measurement value automatically when it judges that the value is stable. Press the MEAS key when the instantaneous value is displayed in the instantaneous measurement, and the HOLD indicator blinks. Then, when the value is stable, the HOLD indicator lights up and the value is held. To cancel the HOLD state or to stop the stability judgment (HOLD blinking), press the MEAS key again.

Stabilization Criteria

pH measurement:

pH change in 10 seconds: within ±0.017 pH (converted to 25°C)

ORP and ion measurements:

Potential change in 10 seconds: within ±1 mV


Conductivity, resistivity, salinity, and TDS measurements: Display change in 10 seconds: within ±3 digit

Temperature measurement:

Temperature change in 10 seconds: within ±2°C

3.4.2 Measurement Items Changing

When the MODE key is pressed, measurement item can be changed in the channel where ACT is lit.




HORIBA

3.4.3 Auto Data Memory

The instrument can perform data memory automatically at a regular interval. Refer to "3.2.3 Auto Data Memory Setting" P.20 for the setting operations.

If the auto power off setting is ON, it is temporarily OFF while this function is active.

 After the interval setting is completed, press the MEAS key to return to the measurement screen.
 "Auto LOG" is displayed on the screen.



2. Press the ENTER key.

"Auto LOGGING" is displayed on the screen, and the auto data memory function starts.



3. To stop the auto data memory, press the ENTER key.

NOTE

- Do not power ON/OFF while the auto data memory function works.
 - It may damage the reliability of the memory data.
- Only the POWER, ENTER, and the ▲ or ▼ keys are possible while the auto data memory function works.
- The actual date of memory may have the difference of ±1 second from the time calculated based on the interval setting.
- When the number of memory data exceeds 999, data memory is stopped and ERROR No.0010 is displayed.


3.4.4 Measurement Value Printing

You can print out measurement values or the HOLD value.

When the automatic printout setting is ON, automatic printout is executed at the following occasions.


- ·You perform data memory during measurement.
- •Calibration or check is completed in the calibration mode.
- **1.** Check that the printer settings and connection are correctly. (Refer to "3.1.4 Connecting and Setting the Printer" P.14.)
- Press the ► key in the measurement mode or HOLD state. The currently displayed measurement value is printed out. Printing examples are shown below.

pH measurement mode



- HOLD condition At HOLD: AUTO Instant value: INST
- ② Temperature compensation setting At Manual mode: MTC At Auto mode: ATC
- ③ Sample ID
- (4) Instrument model
- 5 Serial number
- 6 Electrode status
- ⑦ Asymmetry potential
- (8) Sensitivity

mV measurement mode



- mV status Absolute mV: ABS Relative mV: REL
- ② HOLD condition At HOLD: AUTO Instant value: INST
- ③ Temperature compensation setting At Manual mode: MTC At Auto mode: ATC
- (4) Sample ID
- (5) Instrument model
- 6 Serial number

• ION measurement mode

Date :20	011/03/10	
Time :18	3:39	
ION :17	.7mg/L	_
HOLD :/	AUTO	-1
Temperature	:	
25.0	°C MTC	-2
Sample :00	07	-3
Inst. model :F	-74BW	-4
Inst. SN :12	23456789	-3 -4 -5 -6
Elect. status :	ок —	-6
Offset :	8.1mV	-7)
Sensitivity		
1.55ug/L -1.8	37 g/L	
:100	.0%	<u>_</u> @
1.87 g/L - 20	00 g/L	U
:100	.0%	

- HOLD condition At HOLD: AUTO Instant value: INST
- (2) Temperature compensation setting At Manual mode: MTC At Auto mode: ATC
- 3 Sample ID
- (4) Instrument model
- 5 Serial number
- 6 Electrode status
- ⑦ Asymmetry potential
- (8) Sensitivity

Conductivity measurement mode

- Date :2011/03/21 Time :19:31 COND : 1.90 S/m HOLD :INST-(1) Temperature : 2 25.0°C ATC-3 4 Sample :0000----Inst. model :F-74BW-Inst. SN :123456789 **(5**) CELL 1.014x100 m-1 6 $\overline{(7)}$ Temp. Coef :2.00%/ C
- HOLD condition At HOLD: AUTO Instant value: INST
- ② Temperature compensation setting At Manual mode: MTC At Auto mode: ATC
- 3 Sample ID
- (4) Instrument model
- (5) Serial number
- 6 Cell constant
- ⑦ Temperature coefficient

Salinity measurement mode

Date	:2011/03/10	
Time	:19:22	
SAL	:0.048%	
HOLD	:INST	(1)
Temperat		
	25.0°C MTC	<u> </u>
Sample	:0007	
Inst. mode	el :F-74BW	
Inst. SN	:123456789	
Cal. Coef		6
Temp. Co	ef :2.00%/°C	-+7

- HOLD condition At HOLD: AUTO Instant value: INST
- (2) Temperature compensation setting At Manual mode: MTC At Auto mode: ATC
- 3 Sample ID
- (4) Instrument model
- 5 Serial number
- 6 Calibration coefficient
- ⑦ Temperature coefficient

TDS measurement mode

Date	:2011/03/10	
Time	:19:22	
TDS	: 519mg/L	
HOLD	:INST	
Tempera	iture :	
	25.0°C MTC	2
	:0007	3
Inst. mod	del :F-74BW	
Inst. SN	:123456789	
Temp. Co	oef :2.00%/°C	<u> </u>

- HOLD condition At HOLD: AUTO Instant value: INST
- (2) Temperature compensation setting At Manual mode: MTC At Auto mode: ATC
- 3 Sample ID
- (4) Instrument model
- (5) Serial number
- (6) Temperature coefficient

Resistivity measurement mode

	1
Date :2011/03/10	
Time :19:23	
Resist :	
10.27 ohm m	
HOLD :INST	-(1)
Temperature :	
25.0°C MTC	-2
Sample :0007	-3
Inst. model :F-74BW	-4
Inst. SN :123456789	-5
Temp. Coef :2.00%/ C	-6
-	

- HOLD condition At HOLD: AUTO Instant value: INST
- (2) Temperature compensation setting At Manual mode: MTC At Auto mode: ATC
- 3 Sample ID
- (4) Instrument model
- (5) Serial number
- (6) Temperature coefficient

3.5 Data Operation

3.5.1 Data Memory of Measurement Value

Measurement data can be stored into the instrument memory.


Refer to "3.4.3 Auto Data Memory" P.73 for automatic data memory.

1. Press the DATA key to enter the data operation screen.



2. Press the ENTER key to enter the data memory screen.

The measurement value and "IN" are alternately displayed.



3. Press the ENTER key.

The measurement data is stored into the memory.


3.5.2 Memory Data Display

The memory data can be displayed on the screen.

1. Press the DATA key to enter the data operation mode.




Press the ▲ or ▼ key to enter the memory data display screen.



3. The latest data is displayed.



Press the < or <p>key to switch between date display and DATA No/ID No. display. Press the ▲ or key to search past data.



3.5.3 Memory Data Printing

The memory data can be printed out.

- **1.** Display the memory data referring to "3.5.2 Memory Data Display" P.79.
- **2.** Press the ENTER key.

The displayed measurement data is printed out.

PH measurement mode

Memory Num :001	-1)
Date :2011/03/22	
Time :08:55	
pH : 6.774 HOLD :INST	
HOLD :INST	-2
Temperature :	
25.0°C ATC	-3
Sample :0000	-3 -4
Inst. model :F-74BW	-5
Inst. SN :123456789	-6
Elect. status :OK	-(7)
Memory Num :001	
Date :2011/03/22	
Time :08:55	
COND :0.000mS/m	
HOLD :INST	
Temperature :	
25.0°C MTC	
Sample :0000	
Inst. model :F-74BW	
Inst. SN :123456789	

- (1) Memory number
- ② HOLD condition At HOLD: AUTO Instant value: INST
- Temperature compensation setting At Manual mode: MTC At Auto mode: ATC
- (4) Sample ID
- (5) Instrument model
- 6 Serial number
- ⑦ Electrode status

mV measurement mode

Memory Num :001	-(1)
Date :2011/03/22	
Time :08:55	
mV ABS : 1000.2mV	-2
HOLD :INST	-3
Temperature :	
25.0°C ATC	-(4)
Sample :0000	-4 -5 -6
Inst. model :F-74BW	-6
Inst. SN :123456789	-7
Memory Num :001	
Date :2011/03/22	
Time :08:55	
COND :0.000mS/m	
HOLD :INST	
Temperature :	
25.0°C MTC	
Sample :0000	
Inst. model :F-74BW	
Inst. SN :123456789	

ION measurement mode

Memory Num :001	-(1)
Date :2011/03/22	
Time :08:55	
ION :17.7mg/L	
HOLD :INST	-2
Temperature :	
25.0 C ATC	-3
Sample :0000	-(4)
Inst. model :F-74BW	-3 -4 -5 -6
Inst. SN :123456789	-6
Elect. status :OK	-7
Memory Num :001	
Date :2011/03/22	
Time :08:55	
COND :0.000mS/m	
HOLD :INST	
Temperature :	
25.0°C MTC	
Sample :0000	
Inst. model :F-74BW	
Inst. SN :123456789	

- ① Memory number
- ② mV status Absolute mV: ABS Relative mV: REL
- ③ HOLD condition At HOLD: AUTO Instant value: INST
- Temperature compensation setting At Manual mode: MTC At Auto mode: ATC
- 5 Sample ID
- 6 Instrument model
- ⑦ Serial number
- 1 Memory number
- ② HOLD condition At HOLD: AUTO Instant value: INST
- ③ Temperature compensation setting At Manual mode: MTC At Auto mode: ATC
- (4) Sample ID
- (5) Instrument model
- 6 Serial number
- ⑦ Electrode status


3.5.4 Calibration History Display

• pH calibration history display

1. Press the DATA key to enter the data operation screen.



2. Press the ▲ or ▼ key to enter the calibration history display screen.



3. Press the ENTER key.

The asymmetry potential of the latest calibration data is displayed.



4. Press the **◄** or **▶** key to display the sensitivity and the data of the third point or later.





When the electrode sensitivity is 90% to 93%, "Check" is displayed instead of "Good". When the sensitivity is 85% to 90%, "bAd" is displayed.

Refer to " Electrode status" P.48.

ION calibration history display

1. Press the DATA key to enter the data operation mode.



 Press the ▲ or ▼ key to enter the calibration history display screen.



3. Press the ENTER key.

The asymmetry potential of the latest calibration data is displayed.



4. Press the **◄** or **▶** key to display the sensitivity and the data of the third point or later.




3.5.5 Calibration History Printing

The calibration history can be printed out for each measurement item.

- **1.** Check the printer settings and that the printer is connected correctly. (Refer to "3.1.4 Connecting and Setting the Printer" P.14.)
- **2.** Refer to "3.5.4 Calibration History Display" P.82 to display the calibration history.
- **3.** Press the ENTER key.

The displayed measurement data is printed out.

pH check data (without check data)




- (1) Instrument model
- (2) Serial number
- (3) Electrode state
- (4) Asymmetry potential
- 5 Sensitivity
- 6 Calibration data

PH calibration history (with check data)

- Inst. model :F-74BW $(\mathbf{1})$ Ž 3 Inst. SN :123456789 Elect. state :OK : 0.4mV (4) Offset Sensitivity pH 4.008- 6.865 :100.0% (5) pH 6.865-9.180 : 93.4% Repeatability:Good(0) 0.000 -6 Calibration data Date :2011/03/10 Time :18:56 pH 4.008 : 177.0mV 25.0°C MTC $\overline{(7)}$ pH 6.865 : 7.9mV 25.0°C MTC pH 9.180 : -120.0mV 25.0°C MTC -Inspection before use pH 6.852 : -5.5mV 25.0°C MTC
- ① Instrument model
- Serial number
- ③ Electrode state
- (4) Asymmetry potential
- 5 Sensitivity
- 6 Check value of repeatability
- (7) Calibration data

ION calibration history



- 1 Instrument model
- (2) Serial number
- ③ Electrode state
- (4) Asymmetry potential
- 5 Sensitivity
- 6 Calibration data

3.5.6 pH Check History Display


The pH check history can be printed out.

Calibration Repeatability check

1. Press the DATA key to enter the data operation screen.




2. Press the ▲ or ▼ key to enter the pH check history display screen.



3. Press the ENTER key.

The repeatability of the latest check data is displayed.



3.5.7 pH Check Data Printing

The pH check history can be printed out.

- **1.** Check that the printer settings and that the is connection are correctly. (Refer to "3.1.4 Connecting and Setting the Printer" P.14)
- **2.** Display the memory data referring to "3.5.6 pH Check History Display" P.86.
- **3.** Press the ENTER key.

The displayed measurement data is printed out.

Inst. model :F-74BW	-1
Inst. SN :123456789	-2
Calibration data	
Date :2011/03/10	
Time :18:18	
Inspection befoure use	
Repeatability:Good(0)	
0.000	-3
pH 6.866 : 7.8mV	
25.0°C MTC	-4

- ① Instrument model
- (2) Serial number
- ③ Check value of repeatability
- Temperature compensation setting At Manual mode: MTC At Auto mode: ATC.


Chapter 4 Maintenance and Troubleshooting

This chapter explains the daily maintenance operations and error messages.

Daily maintenance is vital in assuring accurate measurement and preventing breakage from occurring. If electrode maintenance is neglected, it may cause various problems and erroneous measurements. This instrument is equipped with a convenient error message function. If an error message appears, be sure to take appropriate countermeasures.

4.1 pH (ORP) Electrode Maintenance

Maintain the electrodes by referring to the following information or to the operation manuals of the electrodes.




As an example, the maintenance operations of pH electrodes (9615-10D) are described below. Maintain ORP electrodes in the same manner.

Maintenance after daily use


After making measurements, wash the electrode using pure water (deionized), wipe off water from the electrode with filter or tissue paper, and store it with its cap placed on.

Long-term storage




When an electrode is not to be used for a long period of time, store it after performing the following procedure. Also, exchange the internal solution every 3 to 6 months, using the steps explained below.


- **1.** Disconnect the electrode from the instrument.
- **2.** Detach the protective cap from the electrode.



3. Open the internal solution filler port, then use a dropper to remove the internal solution.



4. Fill the electrode with new internal solution (model 300), until it nears the filler port.



5. Wash the tip of the electrode well with pure water (deionized) and wipe it with filter or tissue paper.



6. If the inside of the protective cap is dry, wash the inside with pure water (deionized). After shaking out the water, fill the cap with enough pure water (deionized) to soak the sponge.



If the inside the protective cap of an electrode dries out and the electrode is not used for a long period of time, the response speed of the electrode may slow and its sensitivity may decrease.

Washing electrodes

If the tip of the pH electrode is extremely dirty, response speed may become slow and it may cause errors in measurement value.

If the electrode is so dirty that it cannot be cleaned by rinsing with pure water (deionized), wash the electrode using the methods below depending the dirt.

General dirt & oily grime

Wipe the dirt/grime off using cotton gauze absorbed with neutral detergent.



Inorganic grime

Rinse the tip using a hydrochloric acid solution of approximately 1 mol/L or a wash solution (model 220). Be sure not to soak the electrode in strong acid for a long period of time.



4.2 ION Electrode Maintenance

Refer to the electrode operation manuals for the maintenance operations of each electrode.

lon electrodes

	lon type	Slope(*)	Measurement range	Elect- rode model	Compati ble tip model	Internal solution
Potassium K ⁺	+1	+58 mV	0.04 to 39,000 mg/L	6582 -10C	7682	3.33 mol/L NaCl
Calcium Ca ²⁺	+2	+29 mV	0.40 to 40,080 mg/L	6583 -10C	7683	3.33 mol/L KCI (model 300)
Chloride Cl ⁻	-1	–59 mV	0.4 to 35,000 mg/L	6560 -10C	7660	1 mol/L KNO ₃
Fluoride F ⁻	-1	–59 mV	0.02 to 19,000 mg/L	6561 -10C	7661	3.33 mol/L KCI (model 300)
Nitrate NO ₃ -	-1	–55 mV	0.62 to 62,000 mg/L	6581 -10C	7681	1 mol/L KCl
Ammonia NH ₃	+1	–59 mV	0.1 to 1,000 mg/L	5002A -10C		

(*):

Change in the electric potential of the electrode (25°C) when the ion concentration is changed by a factor of 10.



The above electrodes are subject to change without notice.

The maintenance operations of the ion electrodes listed in the above table are describes below.

Before use

Before using an electrode, condition the electrode according to the following table to prepare it for measurement.

lon electrode	Conditioning agent	Time
Cl⁻ ion electrode	100 mg/L CI- standard solution	Approx. 10 min
F ⁻ ion electrode	Deionized water	Approx. 10 min
NO_3^- ion electrode	1 mol/L potassium nitrate solution (100 g/L KNO ₃)	2 h to 3 h
K ⁺ ion electrode	0.1 mol/L potassium chloride solution (7.5 g/L KCI)	Approx. 12 h
Ca ²⁺ ion electrode	Tap water	Approx. 3 h
NH ₃ ammonia electrode	No conditioning	

Short-term storage

Immerse electrode in the following solutions, when they are to be stored for up to one day and then reused.

lon electrode	Storage solution
Cl ⁻ ion electrode	Deionized water
F⁻ ion electrode	
NO ₃ ⁻ ion electrode	1 mol/L potassium nitrate solution (100 g/L KNO_3)
K ⁺ ion electrode	0.1 mol/L potassium chloride solution (7.5 g/L KCI)
Ca ²⁺ ion electrode	Tap water
NH ₃ ammonia electrode	0.01 mol/L ammonium chloride solution

4.2.1 65XX-10C Electrode Maintenance

Refer to the electrode operation manuals for the maintenance of other electrode models.

Long-term storage

- **1.** Remove the tip electrode from the combination electrode and put on the rubber cap.
- Put on the electrode protective cap.(Do not put water in the electrode protective cap and make sure it is dry.)
- **3.** Store both the tip electrode and the combined electrode in dry condition.
- **4.** To reuse the electrodes, start with the operations explained in " Before use" P.93 above.

Daily maintenance


When an electrode is not used for many hours, the internal solution in the electrode outer tube may be contaminated with samples, or diluted with water. Replace the internal solution in the electrode outer tube once a week or once a month in the following steps.

- **1.** Unplug the rubber stopper from the internal solution filler port, and then turn the electrode upside down and remove the internal solution using a dropper.
- **2.** Use a dropper to fill the electrode outer tube with the specified internal solution.

The ideal amount of internal solution discharge is only a tiny bit from the liquid junction. If the amount of liquid flowing out is extremely small, the electric potential of the reference electrode will not stabilize or will be affected by a stirrer. In such cases, perform the following operations to let the standard solution seep out from the liquid junction.

Internal solution discharge operations

- **1.** Remove the protective tube from the combination electrode, and mount the rubber cap on the tip type ion electrode.
- **2.** Unplug the rubber stopper from the standard solution filler port on the top of the electrode.
- **3.** Stand the electrode vertically on a desk or other horizontal surface, with the bottom of the electrode (the side with the rubber cap) facing down. Push the electrode down two or three times to let the internal solution seep out through the liquid junction.



Filling internal solution (inner tube)

The inner tube of the electrode is of an air-tight structure that allows almost no outflow. Replace the internal solution in the inner tube, however, when the internal solution reduces less than half after a long period of use. (The frequency of internal solution filling is once a year under normal use. It depends on the use and storage conditions.)




To fill the internal solution in the inner tube, you have to disassemble the electrode.

Take extra care for handling the electrode.

Requirements: A dropper, model 330 gel internal solution, the internal solution of the electrode.

Filling the electrode: disassembly procedure


- **1.** Remove the protective tube from the combination electrode, and mount the rubber cap on the tip type ion electrode.
- **2.** Unplug the rubber stopper from the internal solution (outer tube) filler port, and remove the internal solution using a dropper.
- **3.** Screw the electrode cap by hand and pull it 5 cm to 10 cm toward the electrode connector side.




- **4.** Hold the electrode cap and push to stand the electrode vertically on a desk or other horizontal surface, and remove the internal body 1 of the electrode.
- **5.** Move the silicon tube of the internal body 1 downward to expose the internal solution filler port (inner tube).
- **6.** Use a dropper to fill the electrode inner tube with the gel internal solution (model 330) until the solution nears the filler port of the inner tube.

Filling the electrode: assembly

- Return the silicon tube to the original position to seal the internal solution filler port (inner tube).
 Make sure that the filler port is completely sealed.
- **2.** If the bodies 1 or 2 or the liquid junction are dirty, wash them with pure water (deionized).
- Insert the body 1 into the body 2.
 Make sure that the O-ring is securely sealed.
- **4.** Return the spring to the top of body 1.
- **5.** Hold the rubber sleeve in place by hand, then twist the electrode protective cap at 90° and fit the rubber sleeve into the electrode protective cap.



- **6.** Adjust the HORIBA logo on the electrode protective cap and the filling port of the outer tube face in the same direction and fit the cap into the body 2.
- **7.** Twist the rubber sleeve at 90° and hold it in place.
- **8.** Use a dropper to fill the electrode with the specified internal solution.
- **9.** Make the internal solution seep out from the liquid junction referring to " Daily maintenance" P.94.
- **10.**Store the electrode referring to "● Short-term storage" P.93.



4.3 Conductivity Electrode Maintenance

Refer to the electrode operation manuals for the maintenance operations of each electrode.

Long-term storage

When an electrode will not be used for a long period of time, store it after performing the following procedures. Also, perform the maintenance of the electrode every 3 to 6 months.

- **1.** Disconnect the electrode from the instrument.
- **2.** Use pure water to wash away any sample solution that may be adhered to the electrode.
- **3.** Wash the inside of the electrode protective cap with pure water, then, after shaking out the water, fill the cap with enough pure water to soak the sponge.
- **4.** Place the electrode protective cap on the electrode.

4.4 Troubleshooting

This instrument is equipped with a simple error-message function to notify the operator that an operation error or problem with the equipment has occurred. Errors or other problems that occur in the measurement mode are announced by an "ERROR No." appearing in the upper lefthand corner of the screen.

4.4.1 Error Message Chart

ERROR No.	Message	Definition
0001	Memory ERROR	Data cannot be read from or written to the internal memory.
0003	Unstable reading ERROR	The electric potential does not stabilize within three minutes.
0004	Zero calibration ERROR	The asymmetry potential of the electrode is outside the range of ±45 mV.
0005	Span calibration ERROR	pH: The electrode sensitivity is not within the range from 85% to 105% for the theoretical sensitivity. ION: the slope is inverse.
0007	Buffer uncertain	The standard solution can not be identified automatically.
0008	Please calibrate	The actual calibration interval exceeds the calibration cycle setting.
0009	Printer communication ERROR	There is a problem with the printer com- munication.
0010	Data memory over	The number of the memory data exceeds the limit of the memory.
0011	Cell constant out of range	COND: The cell constant is out of the automatic calculation range.
000E	Key operation ERROR	Key operation is invalid.

ERROR No.0001 Memory ERROR

Definition

Data cannot be read from or written to the internal memory.

Cause	Countermeasure
The instrument does not start properly due to noise or other at power ON.	Disconnect and reconnect the AC adapter, and press the ON/OFF key 10 seconds later. If the same error occurs again after the countermeasure is performed, contact HORIBA sales representative.
The internal IC is defective.	Contact HORIBA sales representative.

• ERROR No.0003 Unstable reading ERROR

Definition

The electric potential does not stabilize within three minutes.

Cause	Countermeasure
This is caused by the effect of the sample solution The sample solu- tion is pure water or other solution with low conductivity, or the pH concentra- tion or temperature change.	Press the MEAS key again while "HOLD" blinks or lights up, to measure the sample using instantaneous value measurement.
The electrode is dirty.	Wash the electrode.
The electrode is cracked.	Replace the electrode.
The responsive glass membrane of the electrode has been dry for a long time.	Soak the membrane (on the electrode) in pure water (deionized) for 24 hours.

Cause	Countermeasure
The temperature of the sample solution is fluctuating.	Measure after the sample solution tem- perature stabilizes.

ERROR No.0004 Zero calibration ERROR

Definition

The asymmetry potential of the electrode is outside the range of ± 45 mV.

Cause	Countermeasure
The electrode is dirty.	Wash the electrode.
The electrode is cracked.	Replace the electrode.
The standard solu- tion concentration is fluctuating.	Replace the internal solution in the reference electrode.
The electrode is not connected correctly.	Connect the electrode correctly.
Electrode is not sub- merged deeply enough to cover ref- erence junction.	Immerse the electrode in the sample, to at least 3 cm beyond the electrode tip.
There is problem with the standard buffers.	Use new standard solution.

ERROR No.0005 Span calibration ERROR (pH)

Definition

The electrode sensitivity is not within the range from 85% to 105% for the theoretical sensitivity.

Cause	Countermeasure
The electrode is dirty.	Wash the electrode.
The electrode is cracked.	Replace the electrode.
Calibration was not performed correctly.	Redo the calibration correctly.
There is a problem with the standard solution.	Prepare new standard solutions.
The electrode is not connected correctly.	Connect the electrode correctly.
Electrode is not sub- merged deeply enough to cover ref- erence junction.	Immerse the electrode in the sample, to at least 3 cm beyond the electrode tip.

ERROR No.0007 Buffer uncertain

Definition

The standard solution can not be identified automatically. After performing the countermeasures below, perform calibration again from the beginning.

Cause	Countermeasure
There is a problem with the standard solution.	Prepare new standard solution.
There is a problem with the standard solution setting.	Match the regulation type (NIST, USA) of standard solutions between the settings and the actuals.
The responsive membrane is dry or dirty.	Measure after cleaning the responsive membrane and soaking it in pure water for 24 hours.
Reference internal solution is poor or dirty.	Replace the internal solution with new solution.
The responsive membrane is dam- aged or worn out.	Replace the electrode.

• ERROR No.0008 Please calibrate

Definition

The actual calibration interval exceeds the calibration cycle setting.

Cause	Countermeasure
It exceeds setting days of the calibration cycle.	Perform calibration.

ERROR No.0009 Printer communication ERROR

Definition

There is a problem with the printer communication.

Power OFF the instrument and perform the countermeasures below, and then power ON again.

Cause	Countermeasure
There is a problem with the printer unit connection.	Check the printer connection, and con- nect the instrument and printer again.
The printer is defective.	Replace the printer. If the same error occurs again after the countermeasure is performed, contact HORIBA sales representative.

ERROR No.0010 Data memory over

Definition

The number of the memory data exceeds the limit of the memory.

Cause	Countermeasure
	Check the contents of the memory data, and delete unnecessary data.

• ERROR No.0011 Cell constant out of range

Definition

The cell constant is out of the automatic calculation range.

Cause	Countermeasure
COND electrode is at the end of the useful life.	Replace the electrode.
Improper standard solution	Prepare new standard solution.

• ERROR No.000E Key operation ERROR

Definition

Key operation is invalid.

Cause	Countermeasure
Loose connection or breaking of keys	Disconnect and reconnect the AC adapter, and press the ON/OFF key 10 seconds later. If the same error occurs again after the countermeasure is performed, contact HORIBA sales representative.

4.4.2 More Troubleshooting

This section describes the countermeasures for the troubles that are not indicated by ERROR messages.

• The indicated value fluctuates

When there is a problem with the electrode.

Cause	Countermeasure
The responsive membrane is dry or dirty.	Wash the responsive membrane.
The responsive membrane is dam- aged or worn out.	Replace the electrode.
There are air bub- bles on the elec- trode.	Shake the electrode to remove the air bubbles.
There is no stan- dard solution remaining.	Fill the electrode with new standard solu- tion, as noted in the electrode operation manual.
The wrong stan- dard solution is being used.	Use the correct standard solution.

When there is a problem with the instrument

Cause	Countermeasure
There is a motor or other device caus- ing electrical inter- ference.	Measure at a place where no influence from induction is given. Ground all AC-powered equipment
The electrode is not connected cor- rectly.	Connect the electrode correctly.

Cause	Countermeasure
The liquid junction is not immersed in the sample solu- tion.	Immerse the electrode in the sample solu- tion up until the liquid junction.
Caused by the effects of the sam- ple.	Judge by measuring with a stable stan- dard solution. Contact HORIBA sales representative.

When there is a problem in the sample

• The response is slow.

Cause	Countermeasure
Caused by the effects of the sam- ple.	Response time may slow down, depend- ing on the properties of the sample solu- tion. Contact HORIBA sales representative.
Wash the respon- sive membrane.	Wash the responsive membrane.
The responsive membrane is dam- aged or worn out.	Replace the electrode.
There is a problem with the standard solution.	Fill the electrode with new standard solu- tion, as noted in the electrode operation manual.

The indicated value does not change/No response.

Cause	Countermeasure
Key operation is invalid.	Disconnect and reconnect the AC adapter, and press the ON/OFF key 10 seconds later. If the same error occurs again, seek repairs at contact HORIBA sales representative.
The electrode con- nector is not attached correctly.	Turn the power OFF, then turn it back ON again.
The electrode is defective. (The responsive mem- brane is cracked.)	Replace the electrode.
The pH meter is defective.	Contact HORIBA sales representative.


• The details when the AC adapter is disconnected.

Cause	Countermeasure
No lithium battery is set.	Set the lithium battery (CR2032).
The lithium battery runs down.	Replace the lithium battery (CR2032).

NOTE

Refer to " • Date Setting" P.34, and set date and time.

Lithium battery replacement



• The measured value is blinks.

- •The pH value is outside the measurement range (when pH value is displayed).
 - Measurement range: pH0.000 to pH14.000
- •The mV value is outside the display range (when mV value is displayed).
 - Display range: ±1999.9 mV
- •The ion value is outside the display range (when ion value is displayed).
 - Display range: 99.9 or more
- •The conductivity value is outside the display range (when conductivity value is displayed).

Display range:

0.00 to 19.99 (when cell constant is 100 m^{-1})

Cause	Countermeasure
The sample is inappropriate.	Use a sample within the measurement range.
The liquid junction is not immersed in the sample.	Immerse the electrode up to 3 cm from its tip.
The electrode cable is broken.	Replace the electrode.
The instrument is defective.	Perform the check operation mentioned below.
Calibration is unperformed or incorrect.	Perform calibration correctly.

Check operation

 pH, mV, ion, resistivity Short the connecter for the electrode on the instrument between the center pin and outer metal.



If the measurement value stops blinking in this state, the instrument is normal.

② Conductivity, salinity, TDS Disconnect the electrode. If the displayed measurement value is near 0, the instrument is normal.

The temperature display blinks. The temperature display is fixed at 25°C.

The temperature measurement value is outside the measurement range.

Measurement range: -30°C to 130°C

Cause	Countermeasure
The temperature of the sample is out-side display range.	Check the sample temperature. Adjust it within the measurement range.
The internal ther- mistor connection of the electrode is broken or shorted.	Check the resistance of the temperature sensor connector. If it is 50 k Ω or more at room temperature, replace the electrode.
The connection of the temperature connector is failed.	Insert the temperature connector firmly.
The instrument is defective.	Contact HORIBA sales representative.
The settings of the temperature cali- bration is failed.	Initialize the settings referring to " ● Setting Initialization" P.42.

• Measurements are not repeatable.

Cause	Countermeasure
The effects of the sample.	The pH or other properties of the sample solution may have changed over time, making poor repeatability.
The responsive membrane is dry or dirty.	Wash the responsive membrane.
The internal solu- tion of the elec- trode runs out or contaminated.	Replace the internal solution with new one.
The responsive membrane is cracked or worn out.	Replace the electrode.

• Printout is not made even when a printer is connected.

Check the following points.

- · Is the printer turned ON?
- · Is an error occurred in printer?
- · Is the printing paper run out or clogged?
- When the test print is made according to the manual, is it done correctly?

Chapter 5 Specifications

5.1 Specifications

Measuring object

Measuring object	Item	Description	F-74BW
	Measuring principle	Glass electrode	
	Display range	pH–2.000 to pH19.999	
рН	Measuring range	pH0.000 to pH14.000	\checkmark
	Display resolution	0.001 pH	\checkmark
	Repeatability	±0.005 pH ±1 digit	
Temperature	Measuring principle	Thermistor method	
	Measuring range	0.0°C to 100.0°C	
	Display resolution	0.1°C	\checkmark
	Repeatability	±0.1°C ±1 digit	\checkmark
mV	Measuring range	±1999.9 mV	
	Display resolution	0.1 mV	\checkmark
	Repeatability	±0.1 mV ±1 digit	\checkmark

Measuring object	ltem	Description	F-74BW
ION	Measuring principle	Ion electrode	
	Measuring range	0.00 μg/L to 999 g/L (mol/L)	
	Display resolution	Valid numbers 3 digits	
	Repeatability	±0.5% ±1 digit of full scale	\checkmark
	Measuring principle	2 AC bipolar method	
Conductivity	Measuring range	Cell constant 100 m ⁻¹ 0.000 mS/m to 19.99 S/m Cell constant 10 m ⁻¹ 0.0 μ S/m to 1.999 S/m Cell constant 1000 m ⁻¹ 0.00 mS/m to 199.9 S/m	
	Display resolution	0.05% of full scale	
	Repeatability	±0.5% ±1 digit of full scale	
Resistivity	Measuring principle	Conversion from conductivity value	
	Measuring range	Cell constant 100 m ⁻¹ : 0.00 Ω•m to 199.9 kΩ•m Cell constant 10 m ⁻¹ : 0.0 Ω•m to 1.999 MΩ•m Cell constant 1000 m ⁻¹ : 0.000 Ω•m to 19.99 kΩ•m	\checkmark
	Display resolution	0.05% of full scale	\checkmark
	Repeatability	±0.5% ±1 digit of full scale	\checkmark

Measuring object	Item	Description	F-74BW
Salinity	Measuring principle	Conversion from conductivity value	
	Measuring range	0.00 PPT to 80.00 PPT (0.000% to 8.000%)	
	Display resolution	0.01PPT (0.001%)	
TDS	Measuring principle	Conversion from conductivity value	\checkmark
	Measuring range	0.01 mg/L to 100 g/L	
	Display resolution	0.01 mg/L	\checkmark

• Functions

Item	F-74BW
PC USB interface	
Memory number	999
Calibration points	5
Password setting	
Printer connection	

5.2 Default Settings

Item	Description	Default value	
Common settings	Temperature compensation	Automatic temperature compensation (ATC)	
	Manual temperature compensation value	25°C	
	Auto power off	OFF	
	Sample ID	0000	
	Calibration interval	OFF	
	Password	OFF	
	Auto data memory	OFF	
	Auto print	OFF	
	Calibration standard solution	NIST	
рН	Default value before calibration	Asymmetric potential: 0 mV Sensitivity: 100%	
lon	Ion valency	+1	
	Unit	g/L	
Conductivity	Unit	S/m	
	Temperature coefficient	2.00%/°C (ON)	
	Cell constant	1.000 × 100 m ⁻¹	
Salinity	Unit	PPT	

5.3 Options

This section lists spare and optional parts for the instrument. These parts are possible through HORIBA distributors. Place an order specifying their names, and part numbers.

Part name		Part number	Remarks	
AC adapter	AC adapter, Cable (UL, 120 V)	3014031951		
	AC adapter, Cable (EU, 230 V)	3014031952		
Plain paper printer	Printer (USA 120 V)	3014030146	Printer cable	
	Printer (EU, 230 V)	3014030147	sold separately	
	Printer cable	3014030148	1.5 m	
	Roll paper	3014030149	20 rolls/set	
	Ink ribbon	3014030150	5 pcs/set	
Serial cable		3014030151	1 m	
USB cable		3200373941		
Electrode stand (Standard type)		3200382557		
Electrode stand (Long type)		3200382560		
Stand arm		3200373991		
Sensor holder		3200373961		
X-51 Digital Simulator		_	For pH, mV, ION, DO	
X-52 Digital Simulator		_	For COND	

For any question regarding this product, please contact your local agency, or inquire from the Customer Registration website (http://www.horiba.com/register)

HORIBA Ltd.

1st edition: October 2011 CODE: GZ0000260898