

Magnetic Stirrer


20 Litre Stirrer with Ceramic Plate

Features & Benefits

- Stir up to 20L
- Large square, ceramic glass plate
- Resistant to corrosion and easy to clean
- High magnetic field reduces stirrer bar decoupling
- LCD displays target and actual speeds

Technical Specifications:

Product Code	350025
Plate Dimension	184 x 184mm
Plate Material	Glass Ceramic
Motor Type	DC Brushless Motor
Motor Rating Input	18W
Motor Rating Output	10W
Power	50W
Voltage	240V/60Hz
Stirring Positions	1
Maximum Stirring Volume (H2O)	20L
Maximum Magnetic Stirrer Bar	80mm
Speed Range	100 to 1,500rpm
Speed Display	Digital, LCD
Protection Class	IP21
Dimensions	215 x 360 x 112mm
Weight	5.3g
Permissible Ambient Temperature	5 to 40°C
Permissible Relative Humidity	80%



Ordering Information:

Product Code	Description	Pack Size
350025	Stirrer Ceramic Plate 20L	Each

350025 Magnetic Stirrer Ceramic Plate 20L