

Compact Precision Scale

HT-120

120 g × 0.01 g

HT
SERIES



Stackable Carrying Case
(Standard)



**Handy, Economical,
and Precise!**



AND
A&D Weighing
www.andweighing.com.au

Have You Ever Thought, "It's A Hassle To Maintain Scales"?

If you are tired of ensuring that all scales are properly charged, stored, and ready for use at any time in good condition, you may find the HT-120 compact precision scale to be the perfect solution. Let's take a look:

Unparalleled Battery Life of 500 Hours* * With four AA alkaline batteries, backlight off

- ✓ Saves the trouble of replacing and restocking batteries regularly
- ✓ Auto-power off function to avoid wasting battery power
- ✓ Battery indicator that tells you the remaining battery level
- ✓ Optional AC adapter available

Stainless Steel (SUS304) Weighing Pan and Protector Plate


- ✓ Solid surface that stays free from scratches and wear
- ✓ Protects the scale from solvents, etc.
- ✓ Easy to clean and keep hygienic

- Stackable Carrying Case That Economizes Space for Storage
- Security Slot in the Rear of the Scale to Accommodate an Anti-theft Lock



On top of these features the HT-120 provides first-rate performance and functionality, making it ideal for use in not just schools, but jewelry shops, laboratories and factories as well.

- High Resolution of 0.01 g with Excellent Repeatability
- Large LCD Display (15 mm Character Height) with White Backlight
- Ten Weighing Units: gram, carat, momme, ounce, troy ounce, pennyweight, grain, tola, tael (Singapore/HK jewelry/Taiwan), and Newton
- Tare Function That Offsets the Weight of a Container With Just One Key Press



Specifications

Weighing capacity	120 g / 4.233 oz / 600 ct
Minimum weighing value	0.01 g / 0.001 oz / 0.05 ct
Repeatability	0.01 g (standard deviation)
Linearity	±0.02 g
Power supply	Four AA batteries or AC adapter (optional)

Pan size	Ø70 mm
Dimensions	Main body: 136 (W) × 150 (D) × 57.5 (H) mm Carrying case: 160 (W) × 175 (D) × 70 (H) mm
Net weight	Approx. 570 g (including batteries)
Standard accessories	Instruction manual, carrying case, stainless steel pan, four AA manganese batteries

AC Adapter (Optional)



A&D Australasia
ACN 26 007 556 809
Head Office
32 Dew Street
THEBARTON South Australia 5031
Telephone (08) 8301 8100
Facsimile (08) 8352 7409

Victorian Office
39 Bakehouse Road
KENSINGTON Victoria 3031
Telephone (03) 9371 1555
Facsimile (03) 9372 1193

New South Wales Office
Unit 4, 14 Abbott Road
SEVEN HILLS New South Wales 2147
Telephone (02) 9674 5466
Facsimile (02) 9674 2544

sales@andweighing.com.au
www.andweighing.com.au

Your Dealer

AUSTRALIAN SCIENTIFIC PTY LTD
11 McDOUGALL ST, KOTARA NSW 2289
TEL: 1800 021083 • FAX: 02 49562525
sales@austscientific.com.au • www.austscientific.com.au