



**Washable!**

# **TITAN** Mini Series

Trans-Industrial Technology Advanced Node

## **HL-WP**

Dust & Waterproof Digital Scales

**IP65**



**Stainless Steel Body**



**A&D**  
A&D Australasia  
"Excellence in Measurement"



## Easy & simple to use for sanitary applications!

- Washable body IP65 (stainless steel body)
- Small-footprint
- 1/3000 high resolution  
(HL-300WP, HL-3000WP & HL-3000LWP)
- Portable  
(6 AA-type batteries, batteries not included)
- Stability indicator
- Tare function
- Auto power off
- AC adapter (optional)

### Specifications


| HL-WP                  | HL-300WP                                                       | HL-1000WP | HL-3000WP | HL-3000LWP                                                   |
|------------------------|----------------------------------------------------------------|-----------|-----------|--------------------------------------------------------------|
| Capacity               | 300g                                                           | 1000g     | 3000g     | 3000g                                                        |
| Resolution             | 0.1g                                                           | 0.5g      | 1g        | 1g                                                           |
| Display                | 12mm height, 7 segment Liquid Crystal Display                  |           |           |                                                              |
| Display update         | Approx. 4 times per second                                     |           |           |                                                              |
| Pan size               | 128(W) x 128(D)mm<br>5.04(W) x 5.04(D)inches                   |           |           | 174(W) x 137(D)mm<br>6.85(W) x 5.39(D)inches                 |
| Dimensions             | 170(W) x 220(D) x 63(H)mm<br>6.68(W) x 8.64(D) x 2.50(H)inches |           |           | 174(W) x 220(D) x 63(H)mm<br>6.8(W) x 8.6(D) x 2.50(H)inches |
| Battery life (approx.) | Approx. 100 hours with manganese, 200 hours with alkaline      |           |           |                                                              |
| Weight (approx.)       | Approx. 870g                                                   |           |           | Approx. 950g                                                 |
| Power                  | 6 x R6P / LR6 / "AA" size batteries or AC adapter              |           |           |                                                              |

*Specifications subject to change for improvement without notice.*

### Accessories

AX-TB125 AC Adapter AC200~240V 'C'plug  
AX-TB126 AC Adapter AC200~240V 'BF/S/no'plug  
AX-TB196 AC Adapter AC100V  
AX-TB202 AC Adapter AC120V

### Dimensions (mm)



| Model             | Capacity & Resolution |
|-------------------|-----------------------|
| <b>HL-300WP</b>   | <b>300g x 0.1g</b>    |
| <b>HL-1000WP</b>  | <b>1000g x 0.5g</b>   |
| <b>HL-3000WP</b>  | <b>3000g x 1g</b>     |
| <b>HL-3000LWP</b> | <b>3000g x 1g</b>     |



### TITAN Family

**TITAN** Compact Series  
**SK-WP**  
Dust & Waterproof Digital Scales



| Model                     | Capacity & Resolution |
|---------------------------|-----------------------|
| <b>SK-1000WP</b>          | <b>1000g x 0.5g</b>   |
| <b>SK-2000WP</b>          | <b>2000g x 1g</b>     |
| <b>SK-5000WP</b>          | <b>5000g x 2g</b>     |
| <b>SK-5001WP</b>          | <b>5000g x 1g</b>     |
| <b>SK-10KWP</b>           | <b>10kg x 0.005kg</b> |
| <b>SK-20KWP</b>           | <b>20kg x 0.01kg</b>  |
| .....                     |                       |
| Pan size: 232(W)x192(D)mm |                       |

**AND**  
A&D Australasia

ACN 26 007 556 809

Head Office  
32 Dew Street  
THEBARTON South Australia 5031  
Telephone (08) 8301 8100  
Facsimile (08) 8352 7409

Victorian Office  
39 Bakehouse Road  
KENSINGTON Victoria 3031  
Telephone (03) 9371 1555  
Facsimile (03) 9372 1193

New South Wales Office  
Unit 4, 14 Abbott Road  
SEVEN HILLS New South Wales 2147  
Telephone (02) 9674 5466  
Facsimile (02) 9674 2544

sales@andweighing.com.au  
www.andweighing.com.au

Your Dealer

**AUSTRALIAN SCIENTIFIC PTY LTD**

11 McDOUGALL ST, KOTARA NSW 2289

TEL: 1800 021083 • FAX: 02 49562525

sales@austscientific.com.au • www.austscientific.com.au