

Magnetic Stirrer

20 Litre Stirrer with Ceramic Plate

Features & Benefits


- Stir up to 20L
- Large square, ceramic glass plate
- Resistant to corrosion and easy to clean
- High magnetic field reduces stirrer bar decoupling



Technical Specifications:

Product Code	400.100.021
Plate Dimension	184 x 184mm
Plate Material	Glass Ceramic
Motor Type	DC Brushless Motor
Motor Rating Input	18W
Motor Rating Output	10W
Power	50W
Voltage	240V/60Hz
Stirring Positions	1
Maximum Stirring Volume (H2O)	20L
Maximum Magnetic Stirrer Bar	80mm
Speed Range	100 to 1,500rpm
Protection Class	IP21
Dimensions	215 x 360 x 112mm
Weight	5.3g
Permissible Ambient Temperature	5 to 40°C
Permissible Relative Humidity	80%

400.100.021 Magnetic Stirrer Ceramic Large Plate 20L



Ordering Information:

Product Code	Description	Pack Size
400.100.021	Stirrer Ceramic Plate 20L	Each

AUSTRALIAN  SCIENTIFIC

Australian Scientific Pty Ltd

11 McDougall Street

Kotara NSW 2289

Phone: (02) 4956 2299

Email: sales@austscientific.com.au

Web: www.austscientific.com.au