

AU25302 Poly

150 Litre Corrosive Substance Storage Cabinet



AU25302 Poly Polypropylene Corrosive Cabinet

Complies with AS3780

150 Litre Maximum Capacity

External Size: 1000mm H x 1000mm W x 500mm D

Supplied with two shelves as standard.

The cabinets are fully compliant with the AS 3780 Australian Standards.

Fire Retardant Polypropylene construction Corrosive Cabinets stops the effects of rusting when using highly corrosive vapour acids tested to AS 1530.3—1989.

No Internal Metal Parts.

Adjustable and removable shelving within the cabinet.

Each cabinet is provided standard with a liquid tight sealed sump to contain spills or leaks within the cabinet.

Venting position options are included with the cabinet.

Closing of cabinet - spring loaded door hinge with self-closing doors and locking mechanism.

Opening of cabinet - manual handle opening.

Individual key lockable doors.

12 Months Manufacturers Warranty on parts and labour.


Clean with warm soapy water.

Clean up any spills within the cabinet immediately.

After sales service and training available

AU25302 Poly

150 Litre Corrosive Substance Storage Cabinet



Side View

