

# **AU25452 Poly CW**

## **250 Litre Corrosive Substance Storage Cabinet with Centre Wall Divider for Acids & Basis**



The cabinets are fully compliant with AS 3780 Australian Standards.

Fire Retardant Polypropylene construction Corrosive Cabinets stop the effects of rusting when using highly corrosive vapour acids.

Locking Mechanism, complies with the AS 3780 Australian Standard.

No Internal Metal Parts.

Adjustable and removable shelving within the cabinet.

Each cabinet is provided standard with a liquid tight sealed sump to contain spills or leaks within the cabinet.

Venting position options are included with the cabinet.

Closing of cabinet - spring loaded door hinge with manual slide door handle closing locking.

Opening of cabinet - manual slide handle opening with self closing doors.

Individual key lockable doors.

After sales service and training available

### **AU25452 Poly CW**

Polypropylene Corrosive Cabinet

With Centre Wall Divider

Complies with AS3780

250 Litre Maximum Capacity

External Size: 1760mm H x 1100mm W x 450mm D

# AU25452 Poly CW

## 250 Litre Corrosive Substance Storage Cabinet

WARNING & CAPACITY LABELS

SELF CLOSING HINGES

