

AU25748 Poly CW

100 Litre Corrosive Substance Storage Cabinet With Centre Wall for Acids & Bases



AU25748 Poly CW

Polypropylene Corrosive Cabinet

2 x 50Lt Compartments

Complies with AS3780

100 Litre Maximum Capacity

External Size:

720mm H x 935mm W x 400mm D

Internal Size:

540mm H x 900mm W x 376mm D

Supplied with two shelves as standard.

The cabinets are fully compliant with the relevant Australian Standards.

Manufactured with a Polypropylene centre wall for acids and bases.

Fire Retardant Polypropylene construction Corrosive Cabinets stop the effects of rusting when using highly corrosive vapour acids.

Locking Mechanism, comply to the AS 3780 Australian Standard.

Individual key lockable doors.

No Internal Metal Parts.

Adjustable and removable shelving within the cabinet. With an adjustable distance of 160mm between shelves.

Each cabinet is provided standard with a liquid tight sealed sump to contain spills or leaks within the cabinet.

Venting position options are included with the cabinet.

After sales service and training available

Delivery, most cabinets are available ex stock to meet your requirements.

AU25748 Poly CW
With Centre Wall - Acids & Bases
100 Litre Corrosive Substance Storage Cabinet

