

LABORATORY BLENDERS

1 & 1.2 Liter Blenders

COMMON FEATURES:

Two-speed, push-button operation

Base: die-cast, epoxy-coated motor housing

.5 HP motor



SKU	Technical Information	Approvals	Plug
7011G	120V, 60 Hz	cTUVus	B 
8011G	230V, 50 Hz	RoHS	F 
8011EG	230V, 50 Hz	CE, RoHS	F 

- 1.2-liter, heat-resistant glass container with handle and two-piece vinyl and styrene lid

SKU	Technical Information	Approvals	Plug
7011S	120V, 60 Hz	cTUVus	B 
8011S	230V, 50 Hz	RoHS	F 
8011ES	230V, 50 Hz	CE, RoHS	F 

- 1-liter, stainless steel container with handle and two-piece vinyl and styrene lid