

# LABORATORY BLENDERS



## CB15 SERIES

### Features:

- Easy-to-clean, electronic membrane keypad with 3 speeds plus PULSE
- 4-liter, stainless steel container with two handles and two-piece, easy-off, rubberized lid
- Heavy-duty base with powerful, 3.75 peak input HP motor
- One-piece, dishwasher-safe, removable jar pad

SKU	Technical Info	Approvals
CB15E	230V, 50 Hz	CE, RoHS

*For additional models, see page 23.*

*For full list of accessories, see pages 42-45.*