

Catalogue No: AU25712 Poly

50 Litre Corrosive Substance Storage Cabinet



AU25712 Poly
Polypropylene Corrosive Cabinet
Complies with AS3780
50 Litre Maximum Capacity

External Size:
720mm H x 600mm W x 400mm D
Internal
510mm H x 570mm W x 335mm D

Supplied with one shelf as standard.

The cabinets are fully compliant with the relevant Australian Standards.

Fire Retardant Polypropylene construction Corrosive Cabinets stop the effects of rusting when using highly corrosive vapour acids.

Locking Mechanism, comply to the AS 3780 Australian Standard.

No Internal Metal Parts.

Adjustable and removable shelving within the cabinet.

Each cabinet is provided standard with a liquid tight sealed sump to contain spills or leaks within the cabinet.

Venting position options are included with the cabinet.

Closing of cabinet - spring loaded door hinge with manual slide door handle closing locking.

Opening of cabinet - pull out handle opening with self closing doors.

Individual key lockable doors.

After sales service and training available

Poly Corrosive Substance Storage Cabinets have
(1) Poly shelf for accidental spillages.
Designed to withstand fire the cabinet provide a safe environment to store corrosive chemicals.

