

# Progenics® Nitrile UltraSoft® Powder Free Black Gloves

<b>Product:</b>	Progenics Nitrile <i>UltraSoft</i> ® Powder Free Black Gloves
<b>Material:</b>	Nitrile
<b>Type:</b>	Non Sterile
<b>Design:</b>	Ambidextrous
<b>Feature:</b>	Micro Textured Fingers
<b>Cuff:</b>	Beaded
<b>Usage:</b>	Single Use only
<b>Colour:</b>	Black
<b>Packaging:</b>	XS - L 250/inner, 8 inners/carton XL 240/inner, 8 inners/carton

## Dimensions for UltraSoft Black:

Code	Size	Palm Width (mm)	MGM Length (mm)	ASTM D3578-05 Length (mm)
BPN7851	XS	70 - 79	Min. 240	Min. 220
BPN7852	S	80 - 89	Min. 240	Min. 230
BPN7853	M	90 - 99	Min. 240	Min. 230
BPN7854	L	100 - 109	Min. 240	Min. 230
BPN7855	XL	110 - 119	Min. 240	Min. 230

## Thickness:

Location of Thickness Measurement	Single Wall (mm)
Finger	Min. 0.13
Palm	Min. 0.08
Cuff	Min. 0.05

## Physical Properties:

	Before Aging	After Aging
Tensile Strength (MPa)	Min. 14	Min. 14
Elongation (%)	Min. 500	Min. 400

**Powder Level:** Powder Free



ARTG 174353