

# Centrifuge

## Micro Centrifuge Low Speed

### Features & Benefits

- Ideal for routine medical practice in hospitals, clinics and research institutes
- Built-in fixed angle rotor
- Speed switch (rpm/rcf)
- LCD display and sound alert
- Automatic lid-lock
- Overspeed detection
- Automatic internal diagnosis
- Maintenance free DC brushless motor



### Technical Specifications:

Product Code	400.003.020
Speed Range	300 - 5,000 rpm
Speed Increment / Accuracy	10rpm / $\pm 20$ rpm
Maximum RCF	2,350 x G
Rotor Capacity	6 x 15mL/10mL/7mL/ 1.5-5mL
Run Time	30 sec - 99 min, Continuous
Programs Memory	2
Acceleration/Deceleration Speed	1 $\uparrow$ /2 $\downarrow$
Driving Motor	Brushless DC
Power	Single Phase, 110 - 240 V, 50Hz/60Hz
Dimensions (LxWxH)	300 x 240 x 180 mm
Weight	5.2kg

400.003.020 Micro Centrifuge Low Speed  
(tubes sold separately)

### Ordering Information:

Product Code	Description	Pack Size
400.003.020	Micro Centrifuge Low Speed	Each
<b>Accessories</b>		
400.003.025	Adapter for 5mL Tubes	PK24
400.003.026	Adapter for 1.5/2mL Tubes	PK24