
# Quick Start Guide




# MB90/MB120



### MAKING A SIMPLE MEASUREMENT

Use one of the included absorbent glass fiber pads as a sample to become familiar with the instrument and the various displays.




## ASSEMBLY

CONTROLS

# MB90/MB120





To use it, Click on the top right corner of the interface.

🛅 Default Method name Default Drying profile Standard Drying temperature 105°C Switch-off criterior 10:00 min Display result %MC Back

Place sample pan, close housing and tare



Press the setting items to edit them. Please see the instruction manual for more details.

> For the full instruction manual, please scan the QR code or click the link below:

#### MB120 Instruction Manual:

http://dmx.ohaus.com/WorkArea/showcontent.aspx?id= 4294974363

MB90 Instruction Manual:

http://dmx.ohaus.com/WorkArea/showcontent.aspx?id= 4294974528



P/N 30241177F





MB120 Manual

MB90 Manual

MB90/120 Moisture Analyzers Quick Start Guide

Ohaus Corporation www.ohaus.com