

Plasticware

Pipette Pasteur (Transfer)

Features & Benefits

- Manufactured from non-toxic, Low-Density Polyethylene (LDPE)
- Ideal for transferring and dispensing liquids
- Eliminates the risk of cross contamination
- Excellent transparency
- Uniform wall thickness and precise graduations
- Choice of packaging; bulk, sterile in inner packs of 20, sterile individually wrapped



1mL

1mL Pipette Ordering Information:

| Product Code | Packaging | Sterile | Pack Size |
|--------------|----------------------|---------|-----------|
| 650.050.501 | Bulk | No | PK500 |
| 650.050.511 | Inner packs of 20 | Yes | PK500 |
| 650.050.521 | Individually wrapped | Yes | PK500 |

1mL Technical Specifications:

| | |
|-----------------------|-----------|
| Capacity | 5mL |
| Total Length | 160mm |
| Graduated | Up to 1mL |
| Graduation Increments | 0.25mL |
| Drops per 1mL | 25 |
| Stem Diameter | 5mm |
| Bulb Draw | 3.5mL |



3mL

3mL Pipette Ordering Information:

| Product Code | Packaging | Sterile | Pack Size |
|--------------|----------------------|---------|-----------|
| 650.050.503 | Bulk | No | PK500 |
| 650.050.513 | Inner packs of 20 | Yes | PK500 |
| 650.050.523 | Individually wrapped | Yes | PK500 |

3mL Technical Specifications:

| | |
|-----------------------|-----------|
| Capacity | 7mL |
| Total Length | 160mm |
| Graduated | Up to 3mL |
| Graduation Increments | 0.5mL |
| Drops per 1mL | 21 |
| Stem Diameter | 7.8mm |
| Bulb Draw | 3.5mL |



650.050.501/503
Bulk
Non-sterile



650.050.511/513
Inner Packs of 20
Sterile



650.050.521/523
Individually Wrapped
Sterile