

# ***EJ***

# ***Compact Balances***



- *Unique and eye catching colour design.*
- *10 models available choose the one that suits you.*
- *Your choice of multiple weighing units.*
- *Counting and percentage functions included.*
- *Anti theft lock feature for added security.*
- *Stackable body allows convenient storage for multiple units.*
- *Calculates density automatically using the optional density kit.*

**AND**  
**A&D Australasia**  
*Excellence in Measurement*

**A&D's new compact balances, the EJ Series, offer an outstanding list of features, sure to provide you with excellent value for money.**



- Unique and eye catching colour design.
- Hygienic and rustproof stainless steel weighing pan.
- Provides multiple weighing units, including count and percentages.
- Slot in the rear for an anti theft lock ensuring added security.
- Stackable body allows convenient storage of multiple units.
- Calculates density automatically using the optional density kit.
- Simple keyboard function.
- Battery indicator shows remaining charge.
- Configurable mode key allows customisation of balance to suit users needs.
- Backlit display ideal for poorly lit areas.
- 5 x two decimal place models with round pan (Ø110mm).
- 5 x one decimal place units with square pan (140mm x 127mm).
- Underhook weighing (1 decimal place units only).
- 4 x levelling feet to allow easy set up.
- Comparator function (visual only).
- Modular design for density and breeze break options.

MODEL	EJ-120			EJ-410	EJ-610	EJ-1500			EJ-4100	EJ-6100
Weight capacity	120g			410g	610g	1500g			4100g	6100g
Min. display "d"	0.01g					0.1g				
Weighing pan size	110mm Ø					127mm x 140mm				
Weighing units	g, pcs, %, oz, ozt, ct, mom, dwt, GN, N, SG (density), tola*, tael (HG, HJ, T)*									

\* For special versions only. One of them can be selected and installed at the factory.

### Options

EJ-02: USB interface EJ-03: RS232C interface EJ-07: Underhook for EJ-3000 / EJ-4100 / EJ-6100

EJ-08: Underhook for EJ-1500 / EJ-2000 EJ-11: Breeze break EJ-12: Carrying case

EJ-13: Density determination kit for EJ-120 / EJ-200 / EJ-300 / EJ-410 / EJ-610



**A&D**  
A&D Australasia

"Excellence in Measurement"

ACN 007 556 809

Head Office  
32 Dew Street  
THEBARTON South Australia 5031  
Telephone (08) 8301 8100  
Facsimile (08) 8352 7409

Victorian Office  
39 Bakehouse Road  
KENSINGTON Victoria 3031  
Telephone (03) 9372 1522  
Facsimile (03) 9372 1193

Email: [salesinquiry@andaustralasia.com.au](mailto:salesinquiry@andaustralasia.com.au)  
Internet page: <http://www.andaustralasia.com.au>

New South Wales Office  
Unit 4, 14 Abbott Road  
SEVEN HILLS New South Wales 2147  
Telephone (02) 9674 5466  
Facsimile (02) 9674 2544

Your Dealer

**Australian Scientific Pty Ltd**  
**Phone: 02 49562299**  
**Fax: 02 49562525**  
**Free Call: 1800 021083**  
**Email: [sales@austscientific.com.au](mailto:sales@austscientific.com.au)**