



DRAINING RACK WALL MOUNTING #500080

For storing and draining of different shape glassware. Made of heavy gauge wire (approx. 4mm thick) White Nylon coated, to prevent corrosion. Requires screws for installation (not included).

Size: 570mm width x 600mm height, 150mm depth (approximately)

- Pegs
- Top Row: 55mm length
 - 2nd Row: 95mm length
 - 3rd Row: 55mm length, approx. 20mm width
 - Bottom Row: 165mm length, approx. 20mm width

AUSTRALIAN  SCIENTIFIC

Australian Scientific Pty Ltd

11 McDougall Street

Kotara NSW 2289

Phone: (02) 4956 2299

Email: sales@austscientific.com.au

Web: www.austscientific.com.au