


Metallic Electrode (For ORP Measurement)




Model	Operating Temperature Range (°C)	Electrode Material	Internal Solution	Applications
ORP Electrode 9300-10D Waterproof platinum 3-in-1 type  3014046710 Overall length: 150 mm Diameter of probe: 12 mm Connectors: BNC & phono jack	0-60	Pt / Glass	#300 (KCl)	Waterproof; Platinum on the flat tip allows measurement of small volume samples

Dissolved Oxygen (DO) Electrode & Tips




HORIBA Dissolved Oxygen (DO) electrodes are galvanic probes with integrated temperature sensors. With galvanic DO probes, calibration can be performed immediately and in air. The HORIBA DO probes use unique and innovative tips which are replaceable. No need to replace membranes or refill electrolytes.

Two models are available: a Laboratory model (9520) that can be used for BOD measurements, and a Field immersible model (9551) housed in a rugged casing available in 2m and 10m cable configurations. The Laboratory 9520 DO probe is fitted with a rotor as well as an adaptor to facilitate BOD measurements.

Dissolved Oxygen Electrodes

Model	Measurement Range	Response Time	Temperature Range (°C)	Features
9520-10D For laboratories  3014046711 Overall length: 184 mm Diameter of probe: 15 mm Connectors: BNC & phono jack	0-19.99mg/L DO	20 seconds (90% response time at constant temperature)	0-45	Waterproof; It operates with the built-in temperature sensor and replaceable DO tip 7541.
9551-20D; 9551-100D For field immersible type  3014047090; 3014047091 Overall length: 165 mm Diameter of probe: 32 mm Connectors: BNC & phono jack	0-19.99mg/L DO	30 seconds (90% response time at constant temperature)	0-40	Waterproof; It operates with the built-in temperature sensor and replaceable DO tip 5401.
9552-20D; 9552-50D For field immersible type  3200780939; 3200780941 Overall length: 165 mm Diameter of probe: 30 mm Connectors: BNC & phono jack	0-20.00 mg/L DO	30 seconds (90% response time at constant temperature)	0-50	Waterproof; It operates with the built-in temperature sensor and replaceable DO tip 5402.

Dissolved Oxygen Electrode Tips

Model	Description
7541  3014074145 Overall length: 26.5 mm Diameter: 15 mm	Replacement DO tip for 9520-10D
5401  3014072770	Replacement DO tip for 9551-20D and 9551-100D
5402  3200781553	Replacement DO tip for 9552-20D and 9552-50D