

CORROSIVE STORAGE CABINETS

If you need to store chemicals such as Phenol, Nitric Acid, Sulphuric Acid or other substances with highly corrosive vapours, we suggest you look at the Polyethylene or Polypropylene Cabinets.

These Cabinets are manufactured to comply with AS3780 for the safe storage of Corrosive Substances.



Storage Cabinet
CAT.NO. AU25452 Poly
Capacity 250 Lt
3 Shelves



Storage Cabinet
CAT.NO. AU25302 Poly
Capacity 150 Lt
2 Shelves



Storage Cabinet
AU25712 Poly - Capacity 50 Lt
Single Door, 2 shelves



Storage Cabinet
CAT.NO. AU25748B
Capacity 100 Lt
2 Shelves



Image shows the Centre Wall Only...

Available with Centre Wall (except 50Lt model) for storing of Acids & Bases...

Corrosives Storage Cabinets		Approval	Weight Kg	H x W x D mm
White Polypropylene Construction				
AU25712 Poly	50 Litre, 1 Shelf, 1 Door		25.0	750 x 600 x 400
AU25748 Poly	100 Litre, 2 Shelf, 2 Door		40.0	720 x 935 x 400
AU25748 Poly CW	100 Litre, 2 Compart w/centre wall, 2 Shelf, 2 door		42.0	720 x 935 x 400
AU25302 Poly	150 Litre, 2 Shelf, 2 Door		55.0	1000 x 1000 x 500
AU25380 Poly	200 Litre 3 Shelf, 2 Door		75.0	1250 x 1100 x 450
AU25380 Poly CW	200 Litre, 2 Compart w/centre wall, 3 Shelf, 2 Door		85.0	1250 x 1100 x 450
AU25452 Poly	250 Litre, 3 Shelf, 2 Door		105.0	1760 x 1100 x 450
AU25452 Poly CW	250 Litre, 2 Compartment w/ centre wall 3 Shelf 2 Door		115.0	1760 x 1100 x 450

INNOVATIVE SOLUTIONS FOR THE HANDLING AND STORAGE OF DANGEROUS GOODS