Stirrer Magnetic without Heating

10L | SQUARE PLATE



FEATURES + BENEFITS

- Stirs up to 10L
- Large square, ceramic glass plate (184mm x 184mm)
- Resistant to corrosion and easy to clean
- High magnetic field reduces stirrer bar decoupling
- Maintenance-free DC brushless motor

TECHNICAL SPECIFICATIONS

SPECIFICATIONS	DESCRIPTION
PLATE DIMENSIONS	184 x 184mm
PLATE MATERIAL	Glass Ceramic
MAXIMUM STIRRING VOLUME	10L
MAXIMUM MAGNETIC STIRRER BAR LENGTH	80mm
SPEED RANGE	100–1,500rpm
DISPLAY	-
STIRRING POSITIONS	1
PROTECTION CLASS	IP21
MOTOR TYPE	DC brushless
MOTOR RATING INPUT	18W
MOTOR RATING OUTPUT	10W
POWER	30W
VOLTAGE	240V/60Hz
PERMISSIBLE AMBIENT TEMPERATURE	5–40°C
PERMISSIBLE RELATIVE HUMIDITY	80%
DIMENSIONS (WxDxH)	215 x 360 x 112mm
WEIGHT	3.8kg

Supplied with:

- Power cable
- User manual

Ordering Information

PRODUCT CODE	DESCRIPTION	PACK SIZE
400.100.021	Stirrer Magnetic 10L	Each